

FOOO1262959 - I83G02

BATCH 200-005	
QTY 200 - 016	
SUI0000	
I83G02	
FOOO1262959	
Not Rep Document Water Mark	

Records Centre.  
Ref ID: 332650  
Location: SI 05 237 010

## Pharmaceutical Development / Oral solids and warehousing

## PAGE: 1 of 21

**SCOPE OF THE PREPARATION:** Stability studies and clinical trial

[illegible]

[signature]

Approval for use by the Chief of ORAL SOLIDS and  
WAREHOUSING:

[signature]

## Pharmaceutical Development / Oral solids and warehousing

Page: 2 of 21

**Dosage:** 50 mg (as free base)

[illegible]

**NOTE:** \*DOSE FOR 1 CPS TAKEN CONSIDERING 75-85 AS THE THEORETICAL RATIO % FREE BASE / SALT.

Verifier's signature: \_\_\_\_\_ [signature]

Checked by: \_\_\_\_\_ [signature]

**Pharmacia  
& Upjohn**

Pharmaceutical Development / Oral solids and warehousing

Product: SU10398

Lot: I83G02

Page: 3 of 21

Pharmaceutical form: Capsule

Dosage: 50 mg (as free base)

**CLEANING OF THE EQUIPMENT AND ROOMS**

Once the processing has been completed, clean the processing rooms with: 5% PYRONEG AQUEOUS SOLUTION (METHOD 50/CM/019)

Once the processing has been completed, clean the equipment with: 5% PYRONEG AQUEOUS SOLUTION (METHOD 50/CM/019)

**PROCESSING IDENTIFICATION LABELS**

**CONFORMITY VERIFICATION LABELS**

DATE: 06 / 04 / 01

SIGNATURE: [signature]

**LABELS DELIVERED**

No.: 40

DATE: 12 / 05 / 01

SIGNATURE: [signature]

**ADDITIONAL DELIVERED**

No.: \_\_\_\_\_

DATE: / /

SIGNATURE: \_\_\_\_\_

**LABELS USED**

No.: 35

DATE: 10 / 05 / 01

SIGNATURE: [signature]

**DETERIORATED LABELS**

No.: \_\_\_\_\_

DATE: / /

SIGNATURE: \_\_\_\_\_

**LABELS RETURNED**

No.: 5

DATE: 10 / 05 / 01

SIGNATURE: [signature]

(The returned labels are destroyed)

**LABEL MODEL**

**Pharmacia & Upjohn – Oral Solids Section**

Capsule SU10398 50 mg (as free base)

LOT: I83G02

Prep. Date: 04/2001

FORMULA No.: 83HC01

Date: 12 / 04 / 01

Label No. 8 of 8

[signature]

NOTE: \_\_\_\_\_

Date: 12 / 04 / 01

[signature]

Label No. 16 of 16

**Pharmacia & Upjohn – Oral Solids Section**  
Product: Granulated SU10398 50 mg (as free base)  
LOT: I83K01  
Prep. Date: 04/2001  
FORMULA No.: \_\_\_\_\_

Product: SU10398

Lot: I83G02

Page: 4 of 21

Pharmaceutical form: Capsule

Room: 72

Dosage: 50 mg (as free base)

**WEIGHT VERIFICATION OF  
THE RAW MATERIALS**

DATE	OPER. No.	OPERATION DESCRIPTION	PRODUCTION DATA	OPERATOR	VERIFIER
	1	Check the weight of the active principle/s			
	1/1	PRODUCT: ..... LOT: ..... PRACTICAL WEIGHT: ..... g	Lot: ..... Gross: ..... g Tare: ..... g Net: ..... g Scale ID No.: .....		
	1/2	PRODUCT: ..... LOT: ..... PRACTICAL WEIGHT: ..... g [initials] 06/07/2001	Lot: ..... Gross: ..... g Tare: ..... g Net: ..... g Scale ID No.: .....		
	1/3	PRODUCT: ..... LOT: ..... PRACTICAL WEIGHT: ..... g	Lot: ..... Gross: ..... g Tare: ..... g Net: ..... g Scale ID No.: .....		
		PRODUCT: ..... LOT: ..... PRACTICAL WEIGHT: ..... g	Lot: ..... Gross: ..... g Tare: ..... g Net: ..... g Scale ID No.: .....		

Product: SU10398

Lot: I83G02

Page: 5 of 21

Pharmaceutical form: Capsule

Dosage: 50 mg (as free base)

Room: 72

WEIGHT VERIFICATION OF  
THE RAW MATERIALS

DATE	OPER. No.	OPERATION DESCRIPTION	PRODUCTION DATA	OPERATOR	VERIFIER
01 04 12	2	Check the weight of the following raw materials:			
	2/1	PRODUCT: GRANULATE AT 75% P/P OF SU10398	Lot: I83K01		
		LOT: I83K01	Gross: 7080 g		
		PRACTICAL WEIGHT 4280 g	Tare: 2800 g		
			Net: 4280 g		
			Scale ID No.: 50-BL-32		
	2/2	PRODUCT: T/C OPAQUE SWEDISH ORANGE GEL CPS. Flo 3	Lot: AE283		
		LOT: AE283	Gross: 3450 g		
		PRACTICAL WEIGHT 3185 g	Tare: 260 g		
		[initials] 12/4/01	Net: 3190 g		
			Scale ID No.: 50-BL-32		
	2/3	PRODUCT:	Lot:		
			Gross: g		
		LOT:	Tare: g		
		PRACTICAL WEIGHT g	Net: g		
		Scale ID No.:			
2/4	PRODUCT:	Lot:			
		Gross: g			
	LOT:	Tare: g			
	PRACTICAL WEIGHT g	Net: g			
		Scale ID No.:			
2/5	PRODUCT:	Lot:			
		Gross: g			
	LOT:	Tare: g			
	PRACTICAL WEIGHT g	Net: g			
		Scale ID No.:			

Product: SU10398	Lot: 183G02	Room: _____	Page: <u>6</u> of <u>21</u>
Pharmaceutical form: Capsule	Dosage: 50 mg (as a free base)	DISTRIBUTION into CAPSULES	

DATE	OPER. No.	OPERATION DESCRIPTION	PRODUCTION DATA	OPERATOR	VERIFIER
	<u>3</u>	<u>Distribution into capsules</u>			
01	<u>3/1</u>	Verify the conformity of the hard gelatin shells:			
09		Format No.: <u>3</u>			
12		Body: <u>OPAQUE SWEDISH ORANGE</u>			
		Head: <u>OPAQUE SWEDISH ORANGE</u>			
		Printing: <u>/</u>	Conforms: Yes <input checked="" type="checkbox"/> No <input type="checkbox"/>	[initials]	
	<u>3/2</u>	Weigh 100 empty shells to determine the average weight.	Average shell weight: <u>48.2</u> mg ( $\alpha$ )		
	<u>3/3</u>	Make the <u>ZANASIAZS</u>	Capsule sealing machine: <u>ZANASIAZS</u>	[initials]	[initials]
			ID number: <u>30-09-05</u>		
		type capsule sealing machine ready and set it to	Cleaning verification: <u>OK</u>	[initials]	[initials]
		format No. <u>3</u> with No. <u>2</u> dosage	Dispenser No.: <u>2</u>		
		dispensing.	Format No.: <u>3</u>		
	<u>3/4</u>	Work Parameters <u>89.066</u> [initials] 06/04/01			
		Theoretical weight: <u>89.166</u> mg	Distribution weight: <u>137.26</u> mg ( $\beta$ )	[initials]	[initials]
01		Distribution weight: Theoretical + $\alpha$	Top end weight: <u>143.93</u> mg	[initials]	[initials]
04		Weight limit: $\beta$ + ( $\pm$ <u>7.5</u> % of the theoretical)	Bottom end weight: <u>130.60</u> mg		
12	<u>3/5</u>	Hopper level height: <u>IN PROCESS</u> mm	Hopper level: <u>~30</u> mm	[illegible]	
	<u>3/6</u>	Dispenser chamber height: <u>IN PROCESS</u> mm	Dispenser chamber: <u>11.84</u> mm	VALUE TAKEN	
	<u>3/7</u>	Piston pressure Yes <input checked="" type="checkbox"/> No <input type="checkbox"/>	Pressure index: <u>NOT APPLICABLE</u>	ABOVE THE	
	<u>3/8</u>	Teflon coated pistons Yes <input checked="" type="checkbox"/> No <input type="checkbox"/>	Yes <input checked="" type="checkbox"/> No <input type="checkbox"/>	DISPENSOR	
	<u>3/9</u>	Machine speed: <u>3500</u> cps/h	Machine speed: <u>3500</u> cps/h		
			Production speed: <u>3500</u> cps/h		

Edition No.: 7 of 10/05/99  
Substitutes edition No.: 6 of 03/11/97

Checked by: \_\_\_\_\_ [signature]

	12.04.2004	19.05.2004
0720	0.000	0.000
0721	0.000	0.000
0722	0.000	0.000
0723	0.000	0.000
0724	0.000	0.000
0725	0.000	0.000
0726	0.000	0.000
0727	0.000	0.000
0728	0.000	0.000
0729	0.000	0.000
0730	0.000	0.000
0731	0.000	0.000
0732	0.000	0.000
0733	0.000	0.000
0734	0.000	0.000
0735	0.000	0.000
0736	0.000	0.000
0737	0.000	0.000
0738	0.000	0.000
0739	0.000	0.000
0740	0.000	0.000
0741	0.000	0.000
0742	0.000	0.000
0743	0.000	0.000
0744	0.000	0.000
0745	0.000	0.000
0746	0.000	0.000
0747	0.000	0.000
0748	0.000	0.000
0749	0.000	0.000
0750	0.000	0.000
0751	0.000	0.000
0752	0.000	0.000
0753	0.000	0.000
0754	0.000	0.000
0755	0.000	0.000
0756	0.000	0.000
0757	0.000	0.000
0758	0.000	0.000
0759	0.000	0.000
0760	0.000	0.000
0761	0.000	0.000
0762	0.000	0.000
0763	0.000	0.000
0764	0.000	0.000
0765	0.000	0.000
0766	0.000	0.000
0767	0.000	0.000
0768	0.000	0.000
0769	0.000	0.000
0770	0.000	0.000
0771	0.000	0.000
0772	0.000	0.000
0773	0.000	0.000
0774	0.000	0.000
0775	0.000	0.000
0776	0.000	0.000
0777	0.000	0.000
0778	0.000	0.000
0779	0.000	0.000
0780	0.000	0.000
0781	0.000	0.000
0782	0.000	0.000
0783	0.000	0.000
0784	0.000	0.000
0785	0.000	0.000
0786	0.000	0.000
0787	0.000	0.000
0788	0.000	0.000
0789	0.000	0.000
0790	0.000	0.000
0791	0.000	0.000
0792	0.000	0.000
0793	0.000	0.000
0794	0.000	0.000
0795	0.000	0.000
0796	0.000	0.000
0797	0.000	0.000
0798	0.000	0.000
0799	0.000	0.000
0800	0.000	0.000

0701H	+
0551H	+
0521H	+
0501H	+
0541H	+
0571H	+
0561H	+



Capsules  
Serial Lot

12.04.2001	80025
0	8
0011H	0.007
0021H	0.011
0031H	0.013
0041H	0.015
0051H	0.017
0061H	0.019
0071H	0.021
0081H	0.023
0091H	0.025
0101H	0.027
0111H	0.029
0121H	0.031
0131H	0.033
0141H	0.035
0151H	0.037
0161H	0.039
0171H	0.041
0181H	0.043
0191H	0.045
0201H	0.047
0211H	0.049
0221H	0.051
0231H	0.053
0241H	0.055
0251H	0.057
0261H	0.059
0271H	0.061
0281H	0.063
0291H	0.065
0301H	0.067
0311H	0.069
0321H	0.071
0331H	0.073
0341H	0.075
0351H	0.077
0361H	0.079
0371H	0.081
0381H	0.083
0391H	0.085
0401H	0.087
0411H	0.089
0421H	0.091
0431H	0.093
0441H	0.095
0451H	0.097
0461H	0.099
0471H	0.101
0481H	0.103
0491H	0.105
0501H	0.107
0511H	0.109
0521H	0.111
0531H	0.113
0541H	0.115
0551H	0.117
0561H	0.119
0571H	0.121
0581H	0.123
0591H	0.125
0601H	0.127
0611H	0.129
0621H	0.131
0631H	0.133
0641H	0.135
0651H	0.137
0661H	0.139
0671H	0.141
0681H	0.143
0691H	0.145
0701H	0.147
0711H	0.149
0721H	0.151
0731H	0.153
0741H	0.155
0751H	0.157
0761H	0.159
0771H	0.161
0781H	0.163
0791H	0.165
0801H	0.167
0811H	0.169
0821H	0.171
0831H	0.173
0841H	0.175
0851H	0.177
0861H	0.179
0871H	0.181
0881H	0.183
0891H	0.185
0901H	0.187
0911H	0.189
0921H	0.191
0931H	0.193
0941H	0.195
0951H	0.197
0961H	0.199
0971H	0.201
0981H	0.203
0991H	0.205
1001H	0.207

n		100
043:00	+	0.1
044:00	+	0.1
045:00	+	0.1
046:00	+	0.0
047:00	+	0.0
048:00	+	0.0
049:00	+	0.0
050:00	+	0.0
051:00	+	0.0
052:00	+	0.0
053:00	+	0.0
054:00	+	0.0
055:00	+	0.0
056:00	+	0.0
057:00	+	0.0
058:00	+	0.0
059:00	+	0.0
060:00	+	0.0
061:00	+	0.0
062:00	+	0.0
063:00	+	0.0
064:00	+	0.0
065:00	+	0.0
066:00	+	0.0
067:00	+	0.0
068:00	+	0.0
069:00	+	0.0
070:00	+	0.0
071:00	+	0.0
072:00	+	0.0
073:00	+	0.0
074:00	+	0.0
075:00	+	0.0
076:00	+	0.0
077:00	+	0.0
078:00	+	0.0
079:00	+	0.0
080:00	+	0.0
081:00	+	0.0
082:00	+	0.0
083:00	+	0.0
084:00	+	0.0
085:00	+	0.0
086:00	+	0.0
087:00	+	0.0
088:00	+	0.0
089:00	+	0.0
090:00	+	0.0
091:00	+	0.0
092:00	+	0.0
093:00	+	0.0
094:00	+	0.0
095:00	+	0.0
096:00	+	0.0
097:00	+	0.0
098:00	+	0.0
099:00	+	0.0
100:00	+	0.0
101:00	+	0.0
102:00	+	0.0
103:00	+	0.0
104:00	+	0.0
105:00	+	0.0
106:00	+	0.0
107:00	+	0.0
108:00	+	0.0
109:00	+	0.0
110:00	+	0.0
111:00	+	0.0
112:00	+	0.0
113:00	+	0.0
114:00	+	0.0
115:00	+	0.0
116:00	+	0.0
117:00	+	0.0
118:00	+	0.0
119:00	+	0.0
120:00	+	0.0
121:00	+	0.0
122:00	+	0.0
123:00	+	0.0
124:00	+	0.0
125:00	+	0.0
126:00	+	0.0
127:00	+	0.0
128:00	+	0.0
129:00	+	0.0
130:00	+	0.0
131:00	+	0.0
132:00	+	0.0
133:00	+	0.0
134:00	+	0.0
135:00	+	0.0
136:00	+	0.0
137:00	+	0.0
138:00	+	0.0
139:00	+	0.0
140:00	+	0.0
141:00	+	0.0
142:00	+	0.0
143:00	+	0.0
144:00	+	0.0
145:00	+	0.0
146:00	+	0.0
147:00	+	0.0
148:00	+	0.0
149:00	+	0.0
150:00	+	0.0
151:00	+	0.0
152:00	+	0.0
153:00	+	0.0
154:00	+	0.0
155:00	+	0.0
156:00	+	0.0
157:00	+	0.0
158:00	+	0.0
159:00	+	0.0
160:00	+	0.0
161:00	+	0.0
162:00	+	0.0
163:00	+	0.0
164:00	+	0.0
165:00	+	0.0
166:00	+	0.0
167:00	+	0.0
168:00	+	0.0
169:00	+	0.0
170:00	+	0.0
171:00	+	0.0
172:00	+	0.0
173:00	+	0.0
174:00	+	0.0
175:00	+	0.0
176:00	+	0.0
177:00	+	0.0
178:00	+	0.0
179:00	+	0.0
180:00	+	0.0
181:00	+	0.0
182:00	+	0.0
183:00	+	0.0
184:00	+	0.0
185:00	+	0.0
186:00	+	0.0
187:00	+	0.0
188:00	+	0.0
189:00	+	0.0
190:00	+	0.0
191:00	+	0.0
192:00	+	0.0
193:00	+	0.0
194:00	+	0.0
195:00	+	0.0
196:00	+	0.0
197:00	+	0.0
198:00	+	0.0
199:00	+	0.0
200:00	+	0.0

Product: SU10398	Lot: I83G02	Room: 72	Page: 7 of 21
Pharmaceutical form: Capsule	Dosage: 50 mg (as a free base)	DISTRIBUTION into CAPSULES	

DATE	OPER. No.	OPERATION DESCRIPTION	PRODUCTION DATA	OPERATOR	VERIFIER
01 04 12	<u>3/10</u>	During the distribution, guide the produced capsules into a cyclone separator. <i>AND THEREFORE IN DEPOWDERING. (RECORD THE CHARACTERISTICS IN THE NOTES)</i> [initials] 06/04/01	Model: <u>DS71</u> ID number: <u>50/SC/01</u> Cleaning verification: <u>OK</u> Operative parameters: <u>~ 40-42%</u>	[initials]	[initials]
	<u>3/11</u>	<i>FROM THE DEPOWDERING</i> As they come out of the cyclone, collect the capsules in suitable container/s of: <u>DOUBLE PE BAG AND KRAFT BARREL</u>	Container used: <u>PE BAG AND KRAFT BARREL</u>		
01 04 12	<u>4</u>	<b><u>Sampling and controls</u></b>			
	<u>4/1</u>	Monitor the processing so that the process is executed within the set parameters and perform the following controls according to the manner indicated in SOP SF.CF 004 and the indications shown on the corresponding section of the form.		[initials]	[initials]
01 04 12	<u>5</u>	<b><u>Preparations for the sampling of the finished product for controls</u></b>			
	<u>5/1</u>	At the beginning, middle and end of the distribution into capsules, sample (in a manner equally spread throughout) an overall number of capsules equal to <u>about 600</u> units which are necessary for controls on the finished product.	<input checked="" type="checkbox"/>	[initials]	[initials]
01 04 12	<u>6</u>	<b><u>DISTRIBUTION START</u></b>	<u>12.04.01</u> Date: <u>13:30</u> Time: <u>13:30</u> [initials] 12/04/01	[initials]	[initials]

Product: SU10398	Lot: I83G02	Room: 78	Page: 8 of 21
Pharmaceutical form: Capsule	Dosage: 50 mg (as a free base)	DISTRIBUTION into CAPSULES	

DATE	OPER. No.	OPERATION DESCRIPTION	PRODUCTION DATA	OPERATOR	VERIFIER
01 04 12	Z	<u>Controls while in process</u>			
	7/1	Capsule appearance at the beginning of distribution (that there are no signs of rupture or crushing on the body and/or tips)	Capsule appearance at the beginning of distribution <u>CONFORMS</u>	[initials]	[initials]
	7/2	Uniformity of weight/average weight (SOP SF.CI 051)	<input checked="" type="checkbox"/>		
	7/3	Disintegration (SOP SF.CI 015)	<input checked="" type="checkbox"/>		
	-	Uniformity of contents [initials] <input type="checkbox"/> (sample 30 capsules at the beginning – middle – end of distribution and send the samples to <del>SE/</del> Pharmaceutical Controls) <i>12/04/01</i>	Beginning <input type="checkbox"/> Middle <input type="checkbox"/> End <input type="checkbox"/>		
	7/4	Capsule appearance at the end of distribution (that there are no signs of rupture or crushing on the body and/or tips)	Capsule appearance at the end of distribution <u>CONFORMS</u>	[initials]	[initials]

Product: SU10398	Lot: I83G02	Room: 72	Page: 9 of 21
Pharmaceutical form: Capsule	Dosage: 50 mg (as a free base)	DISTRIBUTION into CAPSULES	

**IN PROCESS WEIGHT CONTROLS (SOP SF.CI 051)**

Scale model: <u>SARTORIUS</u>		ID number: <u>50/BL/37</u>												
Frequency	Avg. theoretical weight	Top end weight	Bottom end weight	No. controls per insp.	No. of Operations									
Start/End of processing/day and every 20 MINUTES	137.26	143.93	130.60	20	7/2									
DATE	TIME	SEE NOTE [initials] 14/05/02 SINGLE WEIGHT VALUES										AVG	S.D.	CV%
12/04/01	13:30	135	136	133	131	135	134	135	136	135	134	134.25	1.33	0.99
		135	135	134	134	134	135	131	135	134	134			
	13:30	132	135	137	136	134	135	135	134	136	135			
		MACHINE STOPPED FOR ADHESION PROBLEMS												
18/04/01	14:07											135.7	1.50	1.11
	14:28											135.5	1.80	1.33
	15:07											136.3	1.80	1.32
19/04/01	8:03											136.6	1.90	1.39
	8:26											137.3	1.80	1.31
	8:54											137.8	2.20	1.60
	9:11	SEE ATTACHED PRINTOUTS										136.5	2.10	1.54
	9:56											136.8	1.60	1.17
	10:20											136.5	2.70	1.98
	10:44											134.2	1.90	1.42
	12:49											136.8	2.00	1.46
	13:13	136	138	137	133	134	139	134	135	134	140	136.7	2.41	1.76
		139	138	138	139	133	140	138	139	136	134			
	13:30											138.0	2.40	1.74
	13:55											136.0	1.80	1.32
	14:12											135.7	2.70	1.99
	14:32											134.6	1.80	1.34
	15:10											133.1	1.90	1.43
	15:49											138.1	5.30	3.84%
20/04/01	7:47											137.6	2.30	1.67
	8:47	SEE ATTACHED PRINTOUTS										140.0	1.80	1.29
	8:56											138.6	1.50	1.08
	9:58											140.0	1.50	1.07

OPERATOR'S SIGNATURE: [signature]

VERIFIER'S SIGNATURE: [signature]

Edition No.: 7 of 10/05/99  
Substitutes edition No.: 6 of 03/11/97

Checked by: [signature]

## Pharmaceutical Development / Oral Solids and Warehousing

Product: SU10398	Lot: I83G02	Room: 72	Page: 1 of 9
Pharmaceutical form: Capsule	Dosage: 50 mg (as a free base)	DISTRIBUTION into CAPSULES	

[illegible]

<b>OPERATOR'S SIGNATURE:</b> _____ [signature] _____	<b>VERIFIER'S SIGNATURE:</b> _____ [signature] _____
Edition No.: 7 of 10/05/99 Substitutes edition No.: 6 of 03/11/97	<b>Checked by:</b> _____ [signature] _____

Product: SU10398	Lot: I83G02	Room: 69	Page: 10 of 21
Pharmaceutical form: Capsule		Dosage: 50 mg (as a free base)	<b>DISTRIBUTION into CAPSULES</b>

**IN PROCESS DISINTEGRATION CONTROLS (SOP SF.CI 015)**

EQUIPMENT: <u>SOTAX DT3 DISINTEGRATOR</u>		ID number: <u>50/DS/01</u>						
FREQUENCY	LIMITS	IMMERSION FLUID	No. controls per inspection	No. OPERATIONS				
Start/End of processing/day and every .../.....	≤ 30 min	TDI WATER 37 °C	6	7/3				
DATE	TIME	CONTROLS ON IMMERSION FLUID			SINGLE VALUES			
12-4-01	13:30	Temp: <u>37.2</u> °C Level: <u>CONFORMS</u>	6'00"	6'45"	7'00"	7'30"	8'10"	8'30"
12-4-01	17:00	Temp: <u>37.2</u> °C Level: <u>CONFORMS</u>	6'00"	6'30"	7'05"	7'20"	8'00"	8'15"
18-4-01	14:00	Temp: <u>37.2</u> °C Level: <u>CONFORMS</u>	6'10"	6'35"	6'55"	7'25"	7'55"	8'20"
18-4-01	17:00	Temp: <u>37.2</u> °C Level: <u>CONFORMS</u>	6'05"	6'25"	6'50"	7'20"	7'50"	8'15"
19-4-01	8:30	Temp: <u>37.2</u> °C Level: <u>CONFORMS</u>	6'00"	6'40"	6'55"	7'15"	7'40"	8'00"
19-4-01	17:00	Temp: <u>37.2</u> °C Level: <u>CONFORMS</u>	6'05"	6'30"	6'55"	7'20"	7'45"	8'15"
20-4-01	8:30	Temp: <u>37.2</u> °C Level: <u>CONFORMS</u>	6'00"	6'20"	6'55"	7'10"	7'40"	8'05"
23-4-01 23-4-04	8:30	Temp: _____ °C Level: <u>CONFORMS</u>	6'50"	7'00"	7'30"	7'55"	8'15"	8'55"
[initials] 23-4-01	14:00	Temp: <u>37.3</u> °C Level: <u>CONFORMS</u>	6'15"	6'45"	7'15"	7'30"	7'55"	8'15"
24-4-01	8:30	Temp: <u>37.4</u> °C Level: <u>CONFORMS</u>	6'10"	6'20"	6'50"	7'20"	7'45"	8'15"
24-4-01	12:00	Temp: <u>37.4</u> °C Level: <u>CONFORMS</u>	6'15"	6'45"	7'00"	7'15"	7'50"	8'20"
		Temp: _____ °C Level: _____		[initials] 24/4/01				

OPERATOR'S SIGNATURE: \_\_\_\_\_ [signature] \_\_\_\_\_

VERIFIER'S SIGNATURE: \_\_\_\_\_ [signature] \_\_\_\_\_

Edition No.: 7 of 10/05/99  
Substitutes edition No.: 6 of 03/11/97

Checked by: \_\_\_\_\_ [signature] \_\_\_\_\_

Product: SU10398	Lot: I83G02	Room: 72	Page: 11 of 21
Pharmaceutical form: Capsule	Dosage: 50 mg (as a free base)	DISTRIBUTION into CAPSULES	

DATE	OPER. No.	OPERATION DESCRIPTION	PRODUCTION DATA	OPERATOR	VERIFIER
24/4 01	8	<b>END OF DISTRIBUTION</b>	Date: 24-4-01 Time: 12:00	[initials]	[initials]
24 /04 /01	2	<b>Controls on the finished product</b>		[initials]	[initials]
	2/1	Using the capsules sampled in point .....5/1....., prepare and carry out the following samples:			
24 /04 /01	2/2	No.: 106 for section controls	No.: 106		
	2/3	No.: 50 for chemical controls	No.: 50		
	2/4	No.: 30 for dissolution and possible technological controls by SF/Pharmaceutical Controls	No.: 30		
	2/5	[initials] 06/04/01 No.: 40 g for bacterial loads	No.: 40 g	[initials]	[initials]
	-	No.: for [initials] 06/04/01	No.: [initials] 06/04/01		
	2/6	COMBINE ANY POSSIBLE EXCESS CAPSULES WITH THE BULK			
24 /04 /01	10	<b>Section technological controls</b>			
	10/1	Perform the following section controls and report the data on the appropriate section regarding the FINISHED PRODUCT <input checked="" type="checkbox"/>	Data reported on: FINISHED PRODUCT <input checked="" type="checkbox"/>		
		DURING PROCESSING <input type="checkbox"/>	DURING PROCESSING <input type="checkbox"/>	[initials]	[initials]
	10/2	Uniformity of weight/average weight (SOP SF.CI 051)	<input checked="" type="checkbox"/>		
	10/3	Disintegration (SOP SF.CI 015)	<input checked="" type="checkbox"/>		



Product: SU10398	Lot: I83G02	Room: 72	Page: 12 of 21
Pharmaceutical form: Capsule	Dosage: 50 mg (as a free base)	DISTRIBUTION into CAPSULES	

DATE	OPER. No.	OPERATION DESCRIPTION	PRODUCTION DATA	OPERATOR	VERIFIER
24 /04 /01	11	<b>Analytical controls on the finished product</b>			
	11/1	Send the above taken samples for the execution of the following controls and fill in the appropriate section regarding the: IN PROCESS ANALYTIC CONTROLS <input type="checkbox"/>  SENDING FOR FINISHED PRODUCT ANALYSIS <input checked="" type="checkbox"/>	Data reported on: IN PROCESS ANALYTIC CONTROLS <input type="checkbox"/> SENDING FOR FINISHED PRODUCT ANALYSIS <input checked="" type="checkbox"/>	[initials]	[initials]
	11/2	Titer <input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		
	11/3	Correlated substances <input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		
24 /04 /01	-	Uniformity of content <input type="checkbox"/>	<input type="checkbox"/>		
	11/4	Karl Fisher <input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		
	11/5	Uniformity of weight <input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		
	11/6	Dissolution <input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		
	11/7	Bacterial load <input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		
	11/8	Other: <u>HPLC IDENTITY</u> <input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	[initials]	[initials]
27 /04 /01	12	<b>Metal detector control</b>			
	12/1	At the end of the distribution, pass the suitable capsules through the metal detector	Model: <u>LOCK METAL DETECTOR</u> ID number: <u>50/MT/04</u> Cleaning verification: <u>Yes</u> Operative parameters: <u>PRESET</u> <u>PROGRAM + VERIFICATION WITH STD</u> <u>SAMPLES - TEST OK.</u>	[initials]	[initials]
	12/2	Verify the number of capsules discarded at the end of the operation	Discarded capsules: Gross: <u>70</u> g Tare: <u>[initials] 20</u> g Net: <u>[initials] 50</u> g Equal to <u>138</u> (number) capsules as calculated based on the average weight 0 <u>367</u>	[initials]	[initials]
27 /04 /01	12/3	Take care to send the discarded capsules to be destroyed.	<input checked="" type="checkbox"/>		

Product: SU10398	Lot: I83G02	Room: 72	Page: 13 of 21
Pharmaceutical form: Capsule	Dosage: 50 mg (as a free base)	DISTRIBUTION into CAPSULES	

DATE	OPER. No.	OPERATION DESCRIPTION	PRODUCTION DATA	OPERATOR	VERIFIER
27 /04 /01	13	<b>Processing yield controls</b>			
	13/1	At the end of the processing collect the capsules and place them in the following primary packaging:  <u>DOUBLE PE BAGS + KRAFT BARREL</u>	Primary packaging used:  <u>DOUBLE PE BAGS + KRAFT BARREL</u>		
	13/2	Determine the quantity of product obtained in ponderal terms.	<u>Ponderal yield:</u> Gross: <u>6000</u> ..... g Tare: <u>290</u> ..... g Net: <u>5710</u> ..... g (H)	[initials]	[initials]
	13/3	Calculate the numeric quantity of the obtained product:	<u>Numeric yield:</u> <u>5710</u> = No. <u>41985</u> (G) <u>0.136 g</u>		
	13/4	Numeric yield = H / average weight <sup>(*)</sup> (*) Obtained by final controls		[initials]	[initials]
27 /04 /01	13/5	End of processing yield: (G / THEORETICAL <sup>(*)</sup> ) * 100 (*) T from page 1	<u>% Yield:</u> <u>41985</u> = <u>87.37</u> % <u>48054</u>		
	14	Calculate the mix quantity and residual shells and see to:  SENDING THE MIX AND SHELLS TO BE DESTROYED <input checked="" type="checkbox"/> SET ASIDE THE MIX <input type="checkbox"/> NOTE: _____	<u>Residual mix</u> <u>Residual shells</u> Gross: <u>69.185 g.</u> Gross: <u>1.300 kg</u> Tare: <u>0.14 g.</u> Tare: <u>0.290 kg</u> Net: <u>69.171 g.</u> Net: <u>1.010 kg</u> <u>69.045 g</u> SENT TO BE DESTROYED <input checked="" type="checkbox"/> SET ASIDE <input type="checkbox"/>	[initials]	[initials]

## IN PROCESS ANALYTICAL CONTROLS

**TO SEND TO FINISHED PRODUCT ANALYSIS**

**Edition No.: 7 of 10/05/99**  
**Substitutes edition No.: 6 of 03/11/97**

Checked by: \_\_\_\_\_ [signature]

Checked by: \_\_\_\_\_ [signature] \_\_\_\_\_



7	150
8	1,100
9	2,000
10	3,000
11	4,000
12	5,000
13	6,000
14	7,000
15	8,000
16	9,000
17	10,000
18	11,000
19	12,000
20	13,000
21	14,000
22	15,000
23	16,000
24	17,000
25	18,000
26	19,000
27	20,000
28	21,000
29	22,000
30	23,000
31	24,000
32	25,000
33	26,000
34	27,000
35	28,000
36	29,000
37	30,000
38	31,000
39	32,000
40	33,000
41	34,000
42	35,000
43	36,000
44	37,000
45	38,000
46	39,000
47	40,000
48	41,000
49	42,000
50	43,000
51	44,000
52	45,000
53	46,000
54	47,000
55	48,000
56	49,000
57	50,000
58	51,000
59	52,000
60	53,000
61	54,000
62	55,000
63	56,000
64	57,000
65	58,000
66	59,000
67	60,000
68	61,000
69	62,000
70	63,000
71	64,000
72	65,000
73	66,000
74	67,000
75	68,000
76	69,000
77	70,000
78	71,000
79	72,000
80	73,000
81	74,000
82	75,000
83	76,000
84	77,000
85	78,000
86	79,000
87	80,000
88	81,000
89	82,000
90	83,000
91	84,000
92	85,000
93	86,000
94	87,000
95	88,000
96	89,000
97	90,000
98	91,000
99	92,000
100	93,000

Product: SU10398	Lot: I83G02	Page: <u>16</u> of <u>21</u>
Pharmaceutical form: Capsule	Dosage: 50 mg (as a free base)	

**TECHNOLOGICAL CONTROLS ON THE FINISHED PRODUCT**

DATE	CONTROL	LIMITS/REFERENCES	RESULT	OPERATOR	VERIFIER						
24 /04 /01	<b>AVERAGE WEIGHT</b>  SOP SF.CI 048	[initials] 24-4-01 Theoretical: <u>137.26</u> <del>136</del> mg Minimum: <sup>139</sup> <u>136.6</u> <del>131</del> mg Maximum: <del>143.93</del> <sup>140.6</sup> <del>626</del> mg [initials] 5/6/01	Average: <u>136.1</u> mg S.D.: <u>2.6</u> C.V.%: <u>1.91</u>								
24 /04 /01	<b>DISINTEGRATION</b>  SOP SF.CI 015	Limit: <u>≤30'</u> Immersion fluid: <u>TDI WATER 37°C</u>  Disks: Yes <input type="checkbox"/> No <input checked="" type="checkbox"/>	Disintegrator: <u>SOTAX OT3</u> ID No.: <u>50</u> / <u>DS</u> / <u>137</u> Immersion fluid: <u>TDI WATER 37°C</u>  Temperature: <u>37.4</u> °C Liquid level: <u>CONFORMS</u> Disks: Yes <input type="checkbox"/> No <input checked="" type="checkbox"/> <table border="1"><tr><td>6'15"</td><td>6'45"</td><td>7'10"</td></tr><tr><td>7'55"</td><td>8'20"</td><td>8'45"</td></tr></table>	6'15"	6'45"	7'10"	7'55"	8'20"	8'45"		
6'15"	6'45"	7'10"									
7'55"	8'20"	8'45"									
	<b>LOSS OF WEIGHT</b>  SOP SF.CI 029	Limit: _____ Temperature: _____ °C Time: _____ [initials] 06/04/01	Equipment: _____ ID No.: _____ / _____ / _____ Temperature: _____ °C Time: _____ minutes Loss of weight: _____ %								
	<b>FRIABILITY</b>  SOP SF.CI 025	Quantity for the control: _____	Friabilimeter: _____ ID No.: _____ / _____ / _____ Initial weight: _____ g Final weight: _____ g Friability: _____ %								

Product: SU10398	Lot: I83G02	Page: <u>17</u> of <u>21</u>
Pharmaceutical form: Capsule	Dosage: 50 mg (as a free base)	

DATE	OPER. No.	OPERATION DESCRIPTION	PRODUCTION DATA	OPERATOR	VERIFIER
27/04 /01	15  15/1	<b><u>Sorting the lot</u></b>  Proceed to the sorting of the sample as described in the following section. At the end of processing, collect a number of samples equal to 3% of the numeric yield of the lots at the end of processing, from various points in the bulk. Report the results on the corresponding page.	Quantity sampled: No. <u>1260</u> [illegible] <u>171 g</u> [initials]	[initials]	[initials]

Edition No.: 7 of 10/05/99  
Substitutes edition No.: 6 of 03/11/97

Checked by: \_\_\_\_\_ [signature] \_\_\_\_\_



Product: SU10398	Lot: I83G02	Page: 18 of 21
Pharmaceutical form: Capsule	Dosage: 50 mg (as a free base)	

**SORTING of the SAMPLES from the PRODUCT OBTAINED at the END OF PROCESSING**

PHARMACEUTICAL FORM: CAPSULE

QUANTITY OBTAINED: No. 41985 (A)

QUANTITY to be SORTED: No. 1260 (B) [Equal to 3% of A]

SORTING LIMITS – PRIMARY DEFECTS: NOT MORE than THREE UNITS

APPEARANCE: \_\_\_\_\_

	DATE 27/4/01	DATE	DATE	DATE	DATE	DATE
LIST OF PRIMARY DEFECTS	NO. OF PIECES	NO. OF PIECES	NO. OF PIECES	NO. OF PIECES	NO. OF PIECES	NO. OF PIECES
CAPSULES BROKEN ON THE TIPS	7					
CAPSULES BROKEN ON THE BODY						
BODY IS VISUALIZED ON THE HEAD						
TOTAL	7					
LIST OF SECONDARY DEFECTS	NO. OF PIECES	NO. OF PIECES	NO. OF PIECES	NO. OF PIECES	NO. OF PIECES	NO. OF PIECES
TOTAL						

NOTE: \_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

Operator's signature [signature]	Verifier's signature [signature]	Checked by: [signature]
-------------------------------------	-------------------------------------	----------------------------

**APPEARANCE CONFORMITY**

LOT CONFORMS for APPEARANCE ☐

LOT DOES NOT CONFORM for APPEARANCE go to UNIT SORTING ☒

SECTION CHIEF SIGNATURE: \_\_\_\_\_ [signature]

# Pharmacia & Upjohn

## Pharmaceutical Development / Oral Solids and Warehousing

Product: SU10398		Lot: I83G02		Page: 19 of 21	
Pharmaceutical form: Capsule		Dosage: 50 mg (as a free base)			
DATE	OPER. No.	OPERATION DESCRIPTION	PRODUCTION DATA	OPERATOR	VERIFIER
09/05/01	15/2	If the results of the sampling sorting are outside the set limits, proceed to unit sorting of the lot as described in the attached form.	<input checked="" type="checkbox"/>	[initials]	[initials]
	15/3	At the end of the sorting operation, send the discarded product to be destroyed.	<input checked="" type="checkbox"/>		
24/04/01	16	<b>Counter sampling</b>			
	16/1	Sample .....100..... (number) units and package them in: <u>P.E. BOTTLES</u>	Quantity sampled: No. <u>100</u> <input checked="" type="checkbox"/>	[initials]	[initials]
09/05/01	17	<b>Final lot yield control</b>			
	17/1	Proceed to the quantitative verification of the available product.	Available product Gross: <u>6.130</u> g Tare: <u>0.435</u> g Net: <u>5.695</u> g (U)	[initials]	[initials]
	17/2	Numeric yield = U / average weight <sup>(*)</sup> (*) Taken from the final controls	Numeric yield = <u>41.875 cps</u> (V) [illegible] [initials] 5.6.01		
10/05/01	17/3	% Yield = (V / THEORETICAL <sup>(*)</sup> ) * 100 (*) T of page 1	<u>10-5-01</u> [initials] % Yield: <u>87.66</u> (Z) 87.14	[initials]	[initials]
10/05/01	18	<b>Deposit in the warehouse</b>			
	18/1	Load the finished product and the counter sample into the SF/Warehouse, stocking them as: <u>27 TA</u> [initials] 3/6/01	<input checked="" type="checkbox"/>	[initials]	[initials]

Edition No.: 7 of 10/05/99  
Substitutes edition No.: 6 of 03/11/97

Checked by: \_\_\_\_\_ [signature]

## Pharmaceutical Development / Oral Solids and Warehousing

## UNIT SORTING

Operator's signature: _____[signature]_____	Verifier's signature: _____[signature]_____
<p>Edition No.: 7 of 10/05/99</p> <p>Substitutes edition No.: 6 of 03/11/97</p>	<p>Checked by: _____[signature]_____</p>

Product: SU10398		Lot: I83G02		Page 20 of 24	
Pharmaceutical form: Capsule		Dosage: 50 mg (as a free base)			
DATE	OPERATION No.	NOTES	OPERATOR	VERIFIER	
		[initials] 14/05/01			
14/05/01	7/2 and 7/3	THE GAPS BETWEEN THE INTERMEDIATE AMONG THE CONTROLS IN WEIGHT GREATER THAN THAT PROVIDED FOR OCCURRED WHERE THE PRODUCTION STOPPED DUE TO THE NEED TO CLEAN THE MACHINE (LOAD HOPPER, DISPENSER, BASIN, ETC.) AND CLEAN OUT THE MIX WHICH WAS PARTICULARLY ADHESIVE. AS THESE STOPS WERE FREQUENT AND PROBLEMATIC THE RESTORATION DECIDED TO RECORD IT IN THE FINAL NOTES. [initials] 14 MAY 2001			
14/05/01	7/2	THE OUT OF LIMIT WEIGHT REVEALED IN THE WEIGHT CONTROL OF: 19 APRIL 2001 AT 15:49 WAS "CORRECTED" BY WEIGHING THE UNIT OF THE CAPSULES RELATIVE TO THE CORRESPONDING PRODUCTION PERIOD. IN PROCESS, CONSIDERING THE PROBLEMS OF THE PROCESS, THE CAPSULES OF VARIOUS PRODUCTION BATCHES WERE SEPARATELY COLLECTED (BETWEEN THE TWO WEIGHT CONTROLS) AND COMBINED WITH THE BULK ONLY AFTER A FINAL CONTROL WAS PERFORMED. [initials] 14 MAY 2001			
	17/2	For the calculation of the yield an average weight equal to 136mg was used rather than 136.1 (the operator incorrectly rounded off the average weight) [signature] 5 June 2001			

Product: SU10398

Lot: I83G02

Page 21 of 21

Pharmaceutical form: Capsule

Dosage: 50 mg (as a free base)

**LOT APPROVAL**

**OPERATIVE VERIFICATION of the "ORAL SOLIDS" SECTION**

NOTES: \_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

SIGNATURE: \_\_\_\_\_ [signature] \_\_\_\_\_ DATE: 14/05/01

**CHIEF of "ORAL SOLIDS and WAREHOUSING" APPROVAL**

RESULTS: APPROVED ☒ REJECTED ☐

NOTES: \_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

SIGNATURE: \_\_\_\_\_ [signature] \_\_\_\_\_ DATE: 15/6/01

**USE AUTHORIZATION OF THE CHIEF of "Q.C./PHARMACEUTICAL CONTROLS"**

RESULTS: APPROVED ☒ REJECTED ☐

NOTES: \_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

SIGNATURE: \_\_\_\_\_ [signature] \_\_\_\_\_ DATE: 15/6/01

SF/ORAL SOLIDS

PRODUCT SU10398 CPS 50 mg (AS FREE BASE)

LOT 183K01

PREPARATION DATE 04/01

ATTACHED INDEXES

1. ☐ ACTIVE PRINCIPLE ANALYSIS REPORT
2. ☐ IN PROCESS ANALYTIC CONTROLS REPORT
3. ☐ PROCESS WATER REPORT
4. ☒ ENVIRONMENTAL PARAMETER MONITORING
5. ☒ RAW MATERIALS/PACKAGING MATERIALS REQUESTS
6. ☒ FINISHED PRODUCT ANALYSIS REPORT
7. ☐ BACTERIAL LOAD REPORT
8. ☒ FINISHED PRODUCT DELIVERY FORM
9. ☒ ANALYSES CERTIFICATE
10. ☒ RAW DATA, in process weight controls.
11. ☒ SCHEDULED DEVIATION 14/01: ONLY INTERMEDIATE CLEANING BETWEEN SU  
[illegible] and SU10398
12. ☐ \_\_\_\_\_
13. ☐ \_\_\_\_\_
14. ☐ \_\_\_\_\_

NOTES: \_\_\_\_\_



[initials] LOT 183G02

	24.04.2001	08:40
n		1
001:H	+	0.138 g
n		2
002:H	+	0.140 g
n		3
003:H	+	0.140 g
n		4
004:H	+	0.141 g
n		5
005:H	+	0.141 g
n		6
006:H	+	0.141 g
n		7
007:H	+	0.137 g
n		8
008:H	+	0.139 g
n		9
009:H	+	0.142 g
n		10
010:H	+	0.139 g
n		11
011:H	+	0.141 g
n		12
012:H	+	0.141 g
n		13
013:H	+	0.138 g
n		14
014:H	+	0.138 g
n		15
015:H	+	0.142 g
n		16
016:H	+	0.138 g
n		17
017:H	+	0.140 g
n		18
018:H	+	0.140 g
n		19
019:H	+	0.137 g
n		20
020:H	+	0.138 g

n	20
$\bar{x}$	0.1396 g
s	0.0016 g
std	1.15 %
Ex	2.791 g
min	0.137 g
max	0.142 g
Diff	0.005 g

	24.04.2001	08:40
n		1
001:H	+	0.139 g
n		2
002:H	+	0.138 g
n		3
003:H	+	0.139 g
n		4
004:H	+	0.135 g
n		5
005:H	+	0.137 g
n		6
006:H	+	0.138 g
n		7
007:H	+	0.138 g
n		8
008:H	+	0.138 g
n		9
009:H	+	0.139 g
n		10
010:H	+	0.138 g
n		11
011:H	+	0.135 g
n		12
012:H	+	0.139 g
n		13
013:H	+	0.138 g
n		14
014:H	+	0.135 g
n		15
015:H	+	0.138 g
n		16
016:H	+	0.135 g
n		17
017:H	+	0.138 g
n		18
018:H	+	0.138 g
n		19
019:H	+	0.138 g
n		20
020:H	+	0.138 g

n	20
$\bar{x}$	0.137 g
s	0.001 g
std	1.0 %
Ex	2.75 g
min	0.135 g
max	0.138 g
Diff	0.003 g

LOT 100100

24.04.2001	09:30
n	1
001:H	+ 0.135 g
n	2
002:H	+ 0.137 g
n	3
003:H	+ 0.137 g
n	4
004:H	+ 0.135 g
n	5
005:H	+ 0.134 g
n	6
006:H	+ 0.134 g
n	7
007:H	+ 0.135 g
n	8
008:H	+ 0.135 g
n	9
009:H	+ 0.134 g
n	10
010:H	+ 0.135 g
n	11
011:H	+ 0.135 g
n	12
012:H	+ 0.137 g
n	13
013:H	+ 0.139 g
n	14
014:H	+ 0.135 g
n	15
015:H	+ 0.133 g
n	16
016:H	+ 0.132 g
n	17
017:H	+ 0.133 g
n	18
018:H	+ 0.135 g
n	19
019:H	+ 0.135 g
n	20
020:H	+ 0.133 g

n	20
$\bar{x}$	0.1352 g
s	0.0018 g
srrel	1.33 %
zk	2.704 g
min	0.132 g
max	0.139 g
Diff	0.007 g

24.04.2001	10:00
n	1
001:H	+ 1.364 g
n	2
002:H	+ 0.135 g
n	3
003:H	+ 0.135 g
n	4
004:H	+ 0.135 g
n	5
005:H	+ 0.134 g
n	6
006:H	+ 0.136 g
n	7
007:H	+ 0.138 g
n	8
008:H	+ 0.135 g
n	9
009:H	+ 0.135 g
n	10
010:H	+ 0.137 g
n	11
011:H	+ 0.135 g
n	12
012:H	+ 0.136 g
n	13
013:H	+ 0.136 g
n	14
014:H	+ 0.135 g
n	15
015:H	+ 0.139 g
n	16
016:H	+ 0.134 g
n	17
017:H	+ 0.137 g
n	18
018:H	+ 0.134 g
n	19
019:H	+ 0.137 g
n	20
020:H	+ 0.135 g

[illegible]

n	20
$\bar{x}$	0.1373 g
s	0.0046 g
srrel	3.34 %
zk	3.946 g
min	0.134 g
max	1.364 g
Diff	1.230 g





n  
 2  
 0.134  
 0.002  
 1.4  
 2.7  
 0.1  
 0.1  
 0.0

23.04.2001 09:  
 001:N + 0.138 1  
 002:N + 0.138 2  
 003:N + 0.141 3  
 004:N + 0.137 4  
 005:N + 0.137 5  
 006:N + 0.139 6  
 007:N + 0.137 7  
 008:N + 0.136 8  
 009:N + 0.136 9  
 010:N + 0.138 10  
 011:N + 0.140 11  
 012:N + 0.135 12  
 013:N + 0.138 13  
 014:N + 0.136 14  
 015:N + 0.136 15  
 016:N + 0.136 16  
 017:N + 0.136 17  
 018:N + 0.136 18  
 019:N + 0.136 19  
 020:N + 0.136 20

n  
 20  
 0.1368  
 0.0018  
 1.32  
 2.735  
 0.132  
 0.140  
 0.0039

23.04.2001 08:44  
 001:N + 0.138 1  
 002:N + 0.133 2  
 003:N + 0.137 3  
 004:N + 0.139 4  
 005:N + 0.137 5  
 006:N + 0.137 6  
 007:N + 0.137 7  
 008:N + 0.140 8  
 009:N + 0.137 9  
 010:N + 0.132 10  
 011:N + 0.137 11  
 012:N + 0.136 12  
 013:N + 0.137 13  
 014:N + 0.139 14  
 015:N + 0.139 15  
 016:N + 0.139 16  
 017:N + 0.136 17  
 018:N + 0.136 18  
 019:N + 0.136 19  
 020:N + 0.137 20

LOT 28362  
20-Apr-2001 7:51:50

20.04.2001	07147
001H	+ 0.135 g
002H	+ 0.136 g
003H	+ 0.137 g
004H	+ 0.140 g
005H	+ 0.139 g
006H	+ 0.139 g
007H	+ 0.135 g
008H	+ 0.138 g
009H	+ 0.142 g
010H	+ 0.138 g
011H	+ 0.134 g
012H	+ 0.141 g
013H	+ 0.136 g
014H	+ 0.135 g
015H	+ 0.139 g
016H	+ 0.137 g
017H	+ 0.136 g
018H	+ 0.136 g
019H	+ 0.142 g
020H	+ 0.136 g

n 20  
x 0.1378 g  
s 0.0028 g  
srel 1.67 %  
Ex 2.732 g  
min 0.134 g  
max 0.142 g  
diff 0.008 g

20.04.2001	08147
001H	+ 0.139 g
002H	+ 0.141 g
003H	+ 0.143 g
004H	+ 0.138 g
005H	+ 0.139 g
006H	+ 0.143 g
007H	+ 0.139 g
008H	+ 0.140 g
009H	+ 0.140 g
010H	+ 0.137 g
011H	+ 0.138 g
012H	+ 0.141 g
013H	+ 0.141 g
014H	+ 0.141 g
015H	+ 0.137 g
016H	+ 0.141 g
017H	+ 0.139 g
018H	+ 0.140 g
019H	+ 0.141 g
020H	+ 0.142 g

n 20  
x 0.1400 g  
s 0.0018 g  
srel 1.28 %  
Ex 2.799 g  
min 0.137 g  
max 0.143 g  
diff 0.006 g

17-AUG-2021 2:06 PM

	19.04.2001	05.05.01
000111	+	3.176 $\frac{1}{2}$
000112	+	3.140 $\frac{1}{2}$
005001	+	3.133 $\frac{1}{2}$
004011	+	3.227 $\frac{1}{2}$
005011	+	3.123 $\frac{1}{2}$
006011	+	3.135 $\frac{1}{2}$
007012	+	3.153 $\frac{1}{2}$
008011	+	3.134 $\frac{1}{2}$
207012	+	3.176 $\frac{1}{2}$
010011	+	3.125 $\frac{1}{2}$
011012	+	3.132 $\frac{1}{2}$
012011	+	3.124 $\frac{1}{2}$
013011	+	3.125 $\frac{1}{2}$
014011	+	3.120 $\frac{1}{2}$
015011	+	3.119 $\frac{1}{2}$
016011	+	3.126 $\frac{1}{2}$
017011	+	3.124 $\frac{1}{2}$
018011	+	3.125 $\frac{1}{2}$
019011	+	3.126 $\frac{1}{2}$
020011	+	3.127 $\frac{1}{2}$

m	20
E	0.1364 g
s	0.0015 g
freq	1.59 s
Tr	2.35 s
min	0.133 g
shc	0.140 g
Diff	0.007 g

19.04.2001		19.04.2001
CO12H	-	0.162 ± 8
CO12H	-	0.178 ± 3
CO21H	+	0.177 ± 4
CO21H	+	0.162 ± 3
CO21H	+	0.176 ± 6
CO21H	+	0.172 ± 7
CO21H	+	0.174 ± 3
CO21H	+	0.177 ± 9
CO21H	+	0.176 ± 10
CO21H	+	0.176 ± 12
CO21H	+	0.176 ± 13
CO21H	+	0.176 ± 15
CO21H	+	0.176 ± 18
CO21H	+	0.176 ± 19
CO21H	+	0.176 ± 20
CO21H	+	0.176 ± 21
CO21H	+	0.176 ± 22
CO21H	+	0.176 ± 23
CO21H	+	0.176 ± 24
CO21H	+	0.176 ± 25
CO21H	+	0.176 ± 26
CO21H	+	0.176 ± 27
CO21H	+	0.176 ± 28
CO21H	+	0.176 ± 29
CO21H	+	0.176 ± 30
CO21H	+	0.176 ± 31
CO21H	+	0.176 ± 32
CO21H	+	0.176 ± 33
CO21H	+	0.176 ± 34
CO21H	+	0.176 ± 35
CO21H	+	0.176 ± 36
CO21H	+	0.176 ± 37
CO21H	+	0.176 ± 38
CO21H	+	0.176 ± 39
CO21H	+	0.176 ± 40
CO21H	+	0.176 ± 41
CO21H	+	0.176 ± 42
CO21H	+	0.176 ± 43
CO21H	+	0.176 ± 44
CO21H	+	0.176 ± 45
CO21H	+	0.176 ± 46
CO21H	+	0.176 ± 47
CO21H	+	0.176 ± 48
CO21H	+	0.176 ± 49
CO21H	+	0.176 ± 50
CO21H	+	0.176 ± 51
CO21H	+	0.176 ± 52
CO21H	+	0.176 ± 53
CO21H	+	0.176 ± 54
CO21H	+	0.176 ± 55
CO21H	+	0.176 ± 56
CO21H	+	0.176 ± 57
CO21H	+	0.176 ± 58
CO21H	+	0.176 ± 59
CO21H	+	0.176 ± 60
CO21H	+	0.176 ± 61
CO21H	+	0.176 ± 62
CO21H	+	0.176 ± 63
CO21H	+	0.176 ± 64
CO21H	+	0.176 ± 65
CO21H	+	0.176 ± 66
CO21H	+	0.176 ± 67
CO21H	+	0.176 ± 68
CO21H	+	0.176 ± 69
CO21H	+	0.176 ± 70
CO21H	+	0.176 ± 71
CO21H	+	0.176 ± 72
CO21H	+	0.176 ± 73
CO21H	+	0.176 ± 74
CO21H	+	0.176 ± 75
CO21H	+	0.176 ± 76
CO21H	+	0.176 ± 77
CO21H	+	0.176 ± 78
CO21H	+	0.176 ± 79
CO21H	+	0.176 ± 80
CO21H	+	0.176 ± 81
CO21H	+	0.176 ± 82
CO21H	+	0.176 ± 83
CO21H	+	0.176 ± 84
CO21H	+	0.176 ± 85
CO21H	+	0.176 ± 86
CO21H	+	0.176 ± 87
CO21H	+	0.176 ± 88
CO21H	+	0.176 ± 89
CO21H	+	0.176 ± 90
CO21H	+	0.176 ± 91
CO21H	+	0.176 ± 92
CO21H	+	0.176 ± 93
CO21H	+	0.176 ± 94
CO21H	+	0.176 ± 95
CO21H	+	0.176 ± 96
CO21H	+	0.176 ± 97
CO21H	+	0.176 ± 98
CO21H	+	0.176 ± 99
CO21H	+	0.176 ± 100

in	20
2	b. 1375
8	0.0015
steel	1.31
Zr	2.745
slm	0.170
upr	0.161
DKIF	0.007


19.04.2002	09:15
00:01	0.137 5
00:02	0.137 5
00:03	0.137 5
00:04	0.137 5
00:05	0.137 5
00:06	0.137 5
00:07	0.137 5
00:08	0.137 5
00:09	0.137 5
00:10	0.137 5
00:11	0.137 5
00:12	0.137 5
00:13	0.137 5
00:14	0.137 5
00:15	0.137 5
00:16	0.137 5
00:17	0.137 5
00:18	0.137 5
00:19	0.137 5
00:20	0.137 5

n	29
r	0.5375
b	0.0772
grob	1.60
La	2.756
eln	0.134
max	0.142
2011	2.008

18.04.2001 15:07

n	ART.ID		
002:N	+	0.136 g	1
n			2
003:N	+	0.137 g	3
n			4
004:N	+	0.140 g	5
n			6
005:N	+	0.134 g	7
n			8
006:N	+	0.137 g	9
n			10
007:N	+	0.135 g	11
n			12
008:N	+	0.138 g	13
n			14
009:N	+	0.138 g	15
n			16
010:N	+	0.138 g	17
n			18
011:N	+	0.135 g	19
n			20
012:N	+	0.133 g	
n			
013:N	+	0.138 g	
n			
014:N	+	0.135 g	
n			
015:N	+	0.134 g	
n			
016:N	+	0.134 g	
n			
017:N	+	0.137 g	
n			
018:N	+	0.136 g	
n			
019:N	+	0.139 g	
n			
020:N	+	0.136 g	
n			
021:N	+	0.136 g	

n	21
$\bar{x}$	0.1363 g
s	0.0018 g
srel	1.32 %
Ex	2.863 g


 LOT 283602  
 All. 40  
 08/29/01

18-Apr-2001 14:13:12

18.04.2001		14:07
n		1
001:N	+	0.136 g
n		2
002:N	+	0.135 g
n		3
003:N	+	0.136 g
n		4
004:N	+	0.135 g
n		5
005:N	+	0.135 g
n		6
006:N	+	0.138 g
n		7
007:N	+	0.137 g
n		8
008:N	+	0.133 g
n		9
009:N	+	0.137 g
n		10
010:N	+	0.135 g
n		11
011:N	+	0.138 g
n		12
012:N	+	0.133 g
n		13
013:N	+	0.135 g
n		14
014:N	+	0.138 g
n		15
015:N	+	0.137 g
n		16
016:N	+	0.135 g
n		17
017:N	+	0.136 g
n		18
018:N	+	0.134 g
n		19
019:N	+	0.134 g
n		20
020:N	+	0.136 g
n		20
x		0.1357 g
s		0.0015 g
srel		1.11 %
Ex		2.713 g
min		0.133 g
max		0.138 g
Diff		0.005 g




LOT 283602

---

18.04.2001	14:28
n	1
001:N	+ 0.136 g
n	2
002:N	+ 0.134 g
n	3
003:N	+ 0.135 g
n	4
004:N	+ 0.135 g
n	5
005:N	+ 0.139 g
n	6
006:N	+ 0.135 g
n	7
007:N	+ 0.135 g
n	8
008:N	+ 0.136 g
n	9
009:N	+ 0.139 g
n	10
010:N	+ 0.134 g
n	11
011:N	+ 0.136 g
n	12
012:N	+ 0.134 g
n	13
013:N	+ 0.134 g
n	14
014:N	+ 0.137 g
n	15
015:N	+ 0.135 g
n	16
016:N	+ 0.138 g
n	17
017:N	+ 0.132 g
n	18
018:N	+ 0.134 g
n	19
019:N	+ 0.137 g
n	20
020:N	+ 0.134 g

n	20
$\bar{x}$	0.1355 g
s	0.0018 g
srel	1.33 %
$\Sigma x$	2.709 g
min	0.132 g
max	0.139 g
Diff	0.007 g

---



~~11~~ LOT 183602

19.04.2001	15:4
n	1
001:N	+ 0.142 g
n	2
002:N	+ 0.141 g
n	3
003:N	+ 0.136 g
n	4
004:N	+ 0.140 g
n	5
005:N	+ 0.138 g
n	6
006:N	+ 0.142 g
n	7
007:N	+ 0.140 g
n	8
008:N	+ 0.138 g
n	9
009:N	+ 0.118 g
n	10
010:N	+ 0.136 g
n	11
011:N	+ 0.136 g
n	12
012:N	+ 0.142 g
n	13
013:N	+ 0.140 g
n	14
014:N	+ 0.138 g
n	15
015:N	+ 0.135 g
n	16
016:N	+ 0.139 g
n	17
017:N	+ 0.140 g
n	18
018:N	+ 0.141 g
n	19
019:N	+ 0.143 g
n	20
020:N	+ 0.136 g

n	20
$\bar{x}$	0.1381 g
s	0.0053 g
srel	3.84 %
$\Sigma x$	2.761 g
min	0.118 g
max	0.143 g
Diff	0.025 g

20-Apr-2001



LOT 253 602  
ATTACHMENT 4 [initials] 29 MAY 2002  
NONE OF THE [illegible] POSSIBLE RISKS OF  
CONTAMINATION OF THE [illegible]  
SEE NOTE [illegible] [initials] 29 MAY 2002

29-05-01 20:46 - 1- F0002 HAN-DRAL  
C4 ALLARM SYSTEM SUMMARY

Point/Action/Event Report with following specifications:  
Start Date/Time : 28-04-01 08:00 Stop Date/Time : 02-05-01 17:00  
Time Range : 1 day 1 hour  
Selected Events : Point Events

Point Name	Point Description
1 45C-23.0-PILOCAL0001	UNID.RIFP.JOINTER H2 LIMITE (25.00)
2 45C-23.0-PILOCAL0002	UNID.RIFP.JOINTER H2 LIMITE (25.00)
3 45C-23.0-PILOCAL0003	UNID.RIFP.JOINTER H2 LIMITE (25.00)
4 45C-23.0-PILOCAL0004	UNID.RIFP.JOINTER H2 LIMITE (25.00)
5 45C-23.0-PILOCAL0005	UNID.RIFP.JOINTER H2 LIMITE (25.00)
6 45C-23.0-PILOCAL0006	UNID.RIFP.JOINTER H2 LIMITE (25.00)
7 45C-23.0-PILOCAL0007	UNID.RIFP.JOINTER H2 LIMITE (25.00)
8 45C-23.0-PILOCAL0008	UNID.RIFP.JOINTER H2 LIMITE (25.00)
9 45C-23.0-PILOCAL0009	UNID.RIFP.JOINTER H2 LIMITE (25.00)
10 45C-23.0-PILOCAL0010	UNID.RIFP.JOINTER H2 LIMITE (25.00)
11 45C-23.0-PILOCAL0011	UNID.RIFP.JOINTER H2 LIMITE (25.00)
12 45C-23.0-PILOCAL0012	UNID.RIFP.JOINTER H2 LIMITE (25.00)
13 45C-23.0-PILOCAL0013	UNID.RIFP.JOINTER H2 LIMITE (25.00)
14 45C-23.0-PILOCAL0014	UNID.RIFP.JOINTER H2 LIMITE (25.00)
15 45C-23.0-PILOCAL0015	UNID.RIFP.JOINTER H2 LIMITE (25.00)
16 45C-23.0-PILOCAL0016	UNID.RIFP.JOINTER H2 LIMITE (25.00)
17 45C-23.0-PILOCAL0017	UNID.RIFP.JOINTER H2 LIMITE (25.00)
18 45C-23.0-PILOCAL0018	UNID.RIFP.JOINTER H2 LIMITE (25.00)
19 45C-23.0-PILOCAL0019	UNID.RIFP.JOINTER H2 LIMITE (25.00)
20 45C-23.0-PILOCAL0020	UNID.RIFP.JOINTER H2 LIMITE (25.00)
21 45C-23.0-PILOCAL0021	UNID.RIFP.JOINTER H2 LIMITE (25.00)
22 45C-23.0-PILOCAL0022	UNID.RIFP.JOINTER H2 LIMITE (25.00)
23 45C-23.0-PILOCAL0023	UNID.RIFP.JOINTER H2 LIMITE (25.00)
24 45C-23.0-PILOCAL0024	UNID.RIFP.JOINTER H2 LIMITE (25.00)
25 45C-23.0-PILOCAL0025	UNID.RIFP.JOINTER H2 LIMITE (25.00)
26 45C-23.0-PILOCAL0026	UNID.RIFP.JOINTER H2 LIMITE (25.00)
27 45C-23.0-PILOCAL0027	UNID.RIFP.JOINTER H2 LIMITE (25.00)
28 45C-23.0-PILOCAL0028	UNID.RIFP.JOINTER H2 LIMITE (25.00)
29 45C-23.0-PILOCAL0029	UNID.RIFP.JOINTER H2 LIMITE (25.00)
30 45C-23.0-PILOCAL0030	UNID.RIFP.JOINTER H2 LIMITE (25.00)
31 45C-23.0-PILOCAL0031	UNID.RIFP.JOINTER H2 LIMITE (25.00)
32 45C-23.0-PILOCAL0032	UNID.RIFP.JOINTER H2 LIMITE (25.00)
33 45C-23.0-PILOCAL0033	UNID.RIFP.JOINTER H2 LIMITE (25.00)
34 45C-23.0-PILOCAL0034	UNID.RIFP.JOINTER H2 LIMITE (25.00)
35 45C-23.0-PILOCAL0035	UNID.RIFP.JOINTER H2 LIMITE (25.00)
36 45C-23.0-PILOCAL0036	UNID.RIFP.JOINTER H2 LIMITE (25.00)
37 45C-23.0-PILOCAL0037	UNID.RIFP.JOINTER H2 LIMITE (25.00)
38 45C-23.0-PILOCAL0038	UNID.RIFP.JOINTER H2 LIMITE (25.00)
39 45C-23.0-PILOCAL0039	UNID.RIFP.JOINTER H2 LIMITE (25.00)
40 45C-23.0-PILOCAL0040	UNID.RIFP.JOINTER H2 LIMITE (25.00)
41 45C-23.0-PILOCAL0041	UNID.RIFP.JOINTER H2 LIMITE (25.00)
42 45C-23.0-PILOCAL0042	UNID.RIFP.JOINTER H2 LIMITE (25.00)
43 45C-23.0-PILOCAL0043	UNID.RIFP.JOINTER H2 LIMITE (25.00)
44 45C-23.0-PILOCAL0044	UNID.RIFP.JOINTER H2 LIMITE (25.00)
45 45C-23.0-PILOCAL0045	UNID.RIFP.JOINTER H2 LIMITE (25.00)
46 45C-23.0-PILOCAL0046	UNID.RIFP.JOINTER H2 LIMITE (25.00)
47 45C-23.0-PILOCAL0047	UNID.RIFP.JOINTER H2 LIMITE (25.00)
48 45C-23.0-PILOCAL0048	UNID.RIFP.JOINTER H2 LIMITE (25.00)
49 45C-23.0-PILOCAL0049	UNID.RIFP.JOINTER H2 LIMITE (25.00)
50 45C-23.0-PILOCAL0050	UNID.RIFP.JOINTER H2 LIMITE (25.00)
51 45C-23.0-PILOCAL0051	UNID.RIFP.JOINTER H2 LIMITE (25.00)
52 45C-23.0-PILOCAL0052	UNID.RIFP.JOINTER H2 LIMITE (25.00)
53 45C-23.0-PILOCAL0053	UNID.RIFP.JOINTER H2 LIMITE (25.00)
54 45C-23.0-PILOCAL0054	UNID.RIFP.JOINTER H2 LIMITE (25.00)
55 45C-23.0-PILOCAL0055	UNID.RIFP.JOINTER H2 LIMITE (25.00)
56 45C-23.0-PILOCAL0056	UNID.RIFP.JOINTER H2 LIMITE (25.00)
57 45C-23.0-PILOCAL0057	UNID.RIFP.JOINTER H2 LIMITE (25.00)
58 45C-23.0-PILOCAL0058	UNID.RIFP.JOINTER H2 LIMITE (25.00)
59 45C-23.0-PILOCAL0059	UNID.RIFP.JOINTER H2 LIMITE (25.00)
60 45C-23.0-PILOCAL0060	UNID.RIFP.JOINTER H2 LIMITE (25.00)
61 45C-23.0-PILOCAL0061	UNID.RIFP.JOINTER H2 LIMITE (25.00)
62 45C-23.0-PILOCAL0062	UNID.RIFP.JOINTER H2 LIMITE (25.00)
63 45C-23.0-PILOCAL0063	UNID.RIFP.JOINTER H2 LIMITE (25.00)
64 45C-23.0-PILOCAL0064	UNID.RIFP.JOINTER H2 LIMITE (25.00)
65 45C-23.0-PILOCAL0065	UNID.RIFP.JOINTER H2 LIMITE (25.00)
66 45C-23.0-PILOCAL0066	UNID.RIFP.JOINTER H2 LIMITE (25.00)
67 45C-23.0-PILOCAL0067	UNID.RIFP.JOINTER H2 LIMITE (25.00)
68 45C-23.0-PILOCAL0068	UNID.RIFP.JOINTER H2 LIMITE (25.00)
69 45C-23.0-PILOCAL0069	UNID.RIFP.JOINTER H2 LIMITE (25.00)
70 45C-23.0-PILOCAL0070	UNID.RIFP.JOINTER H2 LIMITE (25.00)
71 45C-23.0-PILOCAL0071	UNID.RIFP.JOINTER H2 LIMITE (25.00)
72 45C-23.0-PILOCAL0072	UNID.RIFP.JOINTER H2 LIMITE (25.00)
73 45C-23.0-PILOCAL0073	UNID.RIFP.JOINTER H2 LIMITE (25.00)
74 45C-23.0-PILOCAL0074	UNID.RIFP.JOINTER H2 LIMITE (25.00)
75 45C-23.0-PILOCAL0075	UNID.RIFP.JOINTER H2 LIMITE (25.00)
76 45C-23.0-PILOCAL0076	UNID.RIFP.JOINTER H2 LIMITE (25.00)
77 45C-23.0-PILOCAL0077	UNID.RIFP.JOINTER H2 LIMITE (25.00)
78 45C-23.0-PILOCAL0078	UNID.RIFP.JOINTER H2 LIMITE (25.00)
79 45C-23.0-PILOCAL0079	UNID.RIFP.JOINTER H2 LIMITE (25.00)
80 45C-23.0-PILOCAL0080	UNID.RIFP.JOINTER H2 LIMITE (25.00)
81 45C-23.0-PILOCAL0081	UNID.RIFP.JOINTER H2 LIMITE (25.00)
82 45C-23.0-PILOCAL0082	UNID.RIFP.JOINTER H2 LIMITE (25.00)
83 45C-23.0-PILOCAL0083	UNID.RIFP.JOINTER H2 LIMITE (25.00)
84 45C-23.0-PILOCAL0084	UNID.RIFP.JOINTER H2 LIMITE (25.00)
85 45C-23.0-PILOCAL0085	UNID.RIFP.JOINTER H2 LIMITE (25.00)
86 45C-23.0-PILOCAL0086	UNID.RIFP.JOINTER H2 LIMITE (25.00)
87 45C-23.0-PILOCAL0087	UNID.RIFP.JOINTER H2 LIMITE (25.00)
88 45C-23.0-PILOCAL0088	UNID.RIFP.JOINTER H2 LIMITE (25.00)
89 45C-23.0-PILOCAL0089	UNID.RIFP.JOINTER H2 LIMITE (25.00)
90 45C-23.0-PILOCAL0090	UNID.RIFP.JOINTER H2 LIMITE (25.00)
91 45C-23.0-PILOCAL0091	UNID.RIFP.JOINTER H2 LIMITE (25.00)
92 45C-23.0-PILOCAL0092	UNID.RIFP.JOINTER H2 LIMITE (25.00)
93 45C-23.0-PILOCAL0093	UNID.RIFP.JOINTER H2 LIMITE (25.00)
94 45C-23.0-PILOCAL0094	UNID.RIFP.JOINTER H2 LIMITE (25.00)
95 45C-23.0-PILOCAL0095	UNID.RIFP.JOINTER H2 LIMITE (25.00)
96 45C-23.0-PILOCAL0096	UNID.RIFP.JOINTER H2 LIMITE (25.00)
97 45C-23.0-PILOCAL0097	UNID.RIFP.JOINTER H2 LIMITE (25.00)
98 45C-23.0-PILOCAL0098	UNID.RIFP.JOINTER H2 LIMITE (25.00)
99 45C-23.0-PILOCAL0099	UNID.RIFP.JOINTER H2 LIMITE (25.00)
100 45C-23.0-PILOCAL0100	UNID.RIFP.JOINTER H2 LIMITE (25.00)

29-05-01 20:46 - 2- F0002 HAN-DRAL  
C4 ALLARM SYSTEM SUMMARY

LOG TIME	KEYNAME	EVENT ALARM TEXT	VALUE	END UNIT
02-05-01 10:15:11	45C-23.0-PILOCAL0001	UNID.RIFP.JOINTER H2 LIMITE (25.00)	25.0	OK
02-05-01 09:15:42	45C-23.0-PILOCAL0007	UNID.RIFP.JOINTER H2 LIMITE (25.00)	-0.49	FLD
02-05-01 09:15:50	45C-23.0-PILOCAL0007	UNID.RIFP.JOINTER H2 LIMITE (25.00)	-0.22	FLD
04-05-01 02:17:56	45C-23.0-PILOCAL0007	UNID.RIFP.JOINTER H2 LIMITE (25.00)	-0.30	FLD
04-05-01 09:04:00	45C-23.0-PILOCAL0007	UNID.RIFP.JOINTER H2 LIMITE (25.00)	-0.38	FLD
04-05-01 10:17:23	45C-23.0-PILOCAL0007	UNID.RIFP.JOINTER H2 LIMITE (25.00)	0.11	FLD
04-05-01 10:15:01	45C-23.0-PILOCAL0007	UNID.RIFP.JOINTER H2 LIMITE (25.00)	-0.01	FLD
04-05-01 15:24:15	45C-23.0-PILOCAL0007	UNID.RIFP.JOINTER H2 LIMITE (25.00)	-1.16	FLD
04-05-01 15:34:16	45C-23.0-PILOCAL0007	UNID.RIFP.JOINTER H2 LIMITE (25.00)	0.08	FLD
04-05-01 15:16:27	45C-23.0-PILOCAL0007	UNID.RIFP.JOINTER H2 LIMITE (25.00)	-0.01	FLD
04-05-01 15:24:17	45C-23.0-PILOCAL0007	UNID.RIFP.JOINTER H2 LIMITE (25.00)	-0.17	FLD
07-05-01 09:02:49	45C-23.0-PILOCAL0007	UNID.RIFP.JOINTER H2 LIMITE (25.00)	-0.49	FLD
07-05-01 09:02:51	45C-23.0-PILOCAL0007	UNID.RIFP.JOINTER H2 LIMITE (25.00)	-0.38	FLD

29-05-01 20:46 - 3- F0002 HAN-DRAL  
C4 ALLARM SYSTEM SUMMARY

LOG TIME	KEYNAME	EVENT ALARM TEXT	VALUE	END UNIT
07-05-01 10:14:10	45C-23.0-PILOCAL0007	UNID.RIFP.JOINTER H2 LIMITE (25.00)	-0.29	FLD
07-05-01 10:14:12	45C-23.0-PILOCAL0007	UNID.RIFP.JOINTER H2 LIMITE (25.00)	-0.20	FLD

A total of 18 records were found for "Historical Activity Inquiry".

END OF REPORT

SCHEDULED DEVIATION REQUEST	
SECTION: Oral Solid Products R&D	No.: 14/01 (as performed by the QA/Quality Systems Section)
DOCUMENT NUMBER AND TITLE: SF.TD 069 Vers. 2: Cleaning of the equipment for the preparation of oral solids pharmaceutical products.	
PRODUCT/MATERIAL/LOT: SU010398 Lot (A) 5975-MTM-0002 (malate salt of SU011248)	ACTIVITY: Production of granulated lot I83K01 Production of capsules lot I83G02
DESCRIPTION OF PROPOSED DEVIATION: Execution of intermediary cleaning of the equipment previously used for the processing of SU011248 (lots I82K01 and I82G02). The cleaning of the equipment by vacuum and the cleaning of the fluid bed granulator is to be carried out with TDI Water. Then we will proceed to sampling the equipment used in the points indicated in the communication of 2 April (see Attachment 1) solely for informative purposes.	
MOTIVATION: The technical rationale for the deviation is provided in the attached documentation: Attachment 2: Memorandum by Sardar Ali (9/02/01) Attachment 3: Communication by David Hahn (14/02/01)	
SIGNATURE/DATE [signature] 06/04/01	
For a process, provide the start and end dates of the process in advance: 17-20 April 2001 09-30 April 2001 [signature] 06 April 2001	
DEVIATION APPROVAL	
SECTION CHIEF: [signature]	QA/QUALITY SYSTEMS: [signature] April 6th, 2001

Author: Paolo Gatti at itnerpo4  
Date: 4/2/01 4:55 PM  
Priority: Normal  
CC: Rosaria Mariani, Luciano Gambini, Paolo DellaVedova, Mauro Olivieri,  
Donata Giudici at ITNERPO1  
TO: Irma Facchetti

Subject: Re[3]: Intermediate cleaning between the manufacturing of SU011248 and SU010 capsules

Hi Everyone,

Luciano and I have defined which points to sample and analyze in the machines used for the processing of SU11248 cleaned with intermediate cleaning prior to working on SU10398.

It has been decided that one point per machine will be sampled, considering that the result with the greatest residue per unit is superficial after the greater cleaning performed prior to the first lot of SU11248 capsules.

In absolute, the following points will be sampled and analyzed (here is the detailed list for Giorgio who will prepare the swabs accordingly):

Zanasi capsule sealer	Hopper base	(OP/05/1P)
Viani Oscillating Granulator	Rear rotor housing	(GS/03/1P)
Glatt 5 Fluid bed dryer	Spy zone	(LF/02/1P)
Diosna Speed Granulator	Crusher	
Pellegrini V Mixer	Bottom	(MS/27/2P)

I spoke with Giorgio and tomorrow he will take the samples and send the swabs to Rita.

The list of sample points will be inserted into the scheduled deviation that will be drawn up to support the "in campaign" processing of the two products which have different instructions (11248 and 10398).

I will meet tomorrow morning with Luciano and Donata about this.

Bye everyone.

Paolo

---

Reply Separator

Subject: Re[2]: Cleaning intermedio fra mfg capsule SU011248 ed SU010  
Author: Irma Facchetti at itnerpo4  
Date: 02/04/01 14.21

Paolo,  
I'm sorry for the lack of understanding about the deviation.

When operation methods different from those described in a SOP are adopted (such as in this case), it is necessary to follow the procedure regarding the deviations.

Bye,  
Irma

---

Reply Separator

Subject: Re: Cleaning intermedio fra mfg capsule SU011248 ed SU010398  
Author: Paolo Gatti at itnerpo4  
Date: 4/2/01 1:09 PM

Irma,

With regard to the scheduled deviation, I would ask that you forward Sugan's memo and Dave's email to me so that I can get things going as soon as possible with Donata. Only one thing is not clear. It is the first time I have heard about the need for scheduled deviation even though it's been at least a month that we've known we would have only done one intermediate cleaning. I have nothing against doing these documents, and I'm absolutely not arguing, but sometimes it would be better if things were defined a little bit in advance.

In my opinion, the same is true for the sampling and analyses. Tomorrow Giorgio will sample the machines in all the points indicated

by the respective cleaning SOP (I think we all agree on this), so that the analyses can be performed. However the execution time is also linked to the availability of Rita's group, as well as to the decision regarding which points are effectively to be analyzed. I repeat my warning (and I think Rita would agree...) that it is not logical to analyze all the points if they are not strictly necessary according to the rationale with which this verification is to be handled, and which I will evaluate first with Paolo Della Vedova to be sure I have properly understood.

Thank you for your quick update after our chat this morning.

Bye,  
Paolo

---

Reply Separator

Subject: Cleaning intermedio fra mfg capsule SU011248 ed SU010398  
Author: Irma Facchetti at itnerpo4  
Date: 02/04/01 12.49

Paolo,

Speaking as QA, I ask that you:

-open a scheduled deviation request and attach the documents detailing the rationale. There is a Memo by Sugan and an E-mail by Dave which will be formalized into a Memo shortly.

With regard to the sampling requested by Shahe:

-within the bounds of the cleaning procedure, for an intermediate type cleaning, no sampling is provided for. The choice of critical points to be examined will be defined with Paolo Della Vedova.

-it would be advisable to carry out these doses within the shortest possible time provided that there is no data or valid rationale which would allow claims to be made regarding the stability of the product under the conservation conditions and time periods that are to be defined.

Best regards  
Irma




Memorandum

To:	Sardar Ali	From:	Peter Giannousis
Dept:		Dept:	PCPD - Analy & Chem. Dev.
Loc./Tel.		Loc./Tel.	B2-2403 ; X3705
Cc:	Arun Koparkar, James Gage, Bhavesh Patel	Date:	09-Feb-01
Subject:	Genealogy of SU010398 lot (A)5975-MTM-0002		

Dear Sardar,

Per your request, the following is a summary of the genealogy of SU010398 lot (A)5975-MTM-0002.

The reaction scheme that was used to prepare this lot of SU010398 from SU011248 is:



In fact SU011248 lot (A2)5953-TJF-0003 was used as starting material to prepare SU010398 lot (A)5975-MTM-0002. SU010398 is the L-malate salt of SU011248, and as such contains about 75% of SU011248 by weight.

The impurities in SU011248 lot (A2)5953-TJF-0003 were higher than those in the previous lots of SU011248 that were tested in GLP toxicological studies. Therefore SU011248 lot (A2)5953-TJF-0003 was qualified for human use by repeating the 2-week tox study. A memo was issued in early January from Toxicology, certifying that there were no significant differences seen in the tox studies with the new lot versus previous lots of SU011248.

The impurities in SU010398 lot (A)5975-MTM-0002 were found to be similar or lower than those in SU011248 lot (A2)5953-TJF-0003. In fact this lot of SU010398 is being used in 3-month GLP tox studies, with results available in May-June 2001.

Based on these facts, it would be expected that there should be no contamination issues in sequential capsule manufacture, as long as the bulk of the SU011248 and the excipients are removed from the equipment. In other words, one would expect that the API impurities would be comparable, and the amount of freebase left in the equipment should be much less than weighing errors of the L-malate salt.

*Peter Giannousis*  
Peter Giannousis

CONFIDENTIAL

Author: David A Hahn at ITNERPO5

Date: 2/14/01 4:39 PM

Normal

TO: bhavesh-patel@sugen.com at SUGEN, chandu-hegde@sugen.com at SUGEN,  
peter-giannousis@sugen.com at SUGEN, sardar-ali@sugen.com at SUGEN

CC: Marco Adami at ITNERPO4, Marina Baldi at ITNERPO4, Irma Facchetti at ITNERPO4,  
Paolo Gatti at ITNERPO4, Rosaria Mariani at ITNERPO4, Mauro Olivieri at ITNERPO4,  
Luciano Gambini at ITNERPO4

-Subject: Re: Sequential capsule manufacturing from free base and L-ma  
----- Message Contents

Sardar,

Here is the general logic that I have in mind. This could be developed in more depth, or in a different way (to the extent allowed by the data). Please let me know what you think.

(1) Solubility and rotating disk dissolution rate data indicate that both the free base and the malate salt have solubility "more than sufficient to prevent solubility from being a limiting factor in the bioavailability," according to Study Report a0089789. Thus, a small amount of one material in the other would not be expected to have any impact on biological performance.

(2) Paolo estimates that after the proposed dry cleaning that the amount of granulation remaining would certainly be less than 10 grams (probably much less). If as much as 10 g remained, this would amount to less than 0.2% of a 5.3 kg granulation batch (using 3.5 kg FBE of API). Given the similar dissolution behaviors, the presence of 0.2% of a granulation of one salt in a granulation of the other would not be expected to influence the biological performance.

(3) Because the process uses wet granulation, the granulated material of which traces would remain on the surfaces of the equipment would likely be representative of the previous granulation, and would likely be incorporated homogeneously into the subsequent granulation. Thus, the presence of a small amount of material from the previous granulation would not be expected to significantly alter the chemical or physical properties of the subsequent granulation.

(4) And I understand that Peter is developing a rationale for safety of the impurity levels based on the genealogical relationship between the batches involved and based on the fact that the qualified impurities levels would allow use of either batch in humans.

Please let me know if you have any comments, questions or concerns.

Ciao,

Dave

Reply Separator

Subject: Sequential capsule manufacturing from free base and L-malate

Author: Sardar Ali <sardar-ali@sugen.com> at SMTP-KZO

Date: 2/12/01 5:37 PM

Dear Dr. Hahn,

Reference to our meeting in Nerviano (during my visit) with Irma, and Rita dated 1/25/01. We had discussed the impact of sequential capsule manufacturing from free base and L-malate salt API's. As we discussed that the equipment will be dry cleaned (removing excipients from the equipment) after completing one batch and before moving to the next batch with different excipient. What we agreed was to get some scientific rationale from you and Peter to assure that the amount of free base traces left in the equipment will be non-detectable. Peter is preparing a summary of genealogy of SU010398 Lot (A) 5975-MTM-002 to justify that impurities in SU010398 are similar or lower than those in SU011248 lot. I will appreciate if we get some scientific rationale from you regarding this what you had agreed to provide us.

I am sorry that I did not get back to you earlier because the manufacturing plan was changed when I returned (Capsule manufacturing from the free base API only) but now it has been changed back to the same what we had discussed With Best Regards

Sardar Ali

QA Product Release Manager

SUGEN, Inc.

230 East Grand Avenue

South San Francisco

Phone: (650) 837-3648

Fax: (650) 837-3326

LOT I83K01 [initials] 29 MAY 2002  
ATTACHMENT 5

## PHARMACEUTICAL DEVELOPMENT

### FINISHED PRODUCT AND PACKAGING MATERIAL REQUEST

<u>SU10398</u>	PRODUCT:	PHARMACEUTICAL FORM: <u>GRANULATE</u>
<u>75% P/P</u>	Dosage:	Lot: <u>I83K01</u>
Approved: <input type="radio"/>	Under analysis: <input checked="" type="radio"/>	
Quantity requested: <u>L280</u> units equal to <u>/</u> g (Average unit weight: <u>)</u>		
Quantity sent: <u>4280*</u> units _____		
Product:	<u>loose x packaged</u>	
Purpose of the request:	<u>stability packaging</u>	
	other: <u>Preparation lot I83G02</u>	

PACKAGING MATERIAL					
MATERIAL	CODE	LOT	QUANTITY REQUESTED	QUANTITY SENT	UNDER ANALYSIS

REQUESTING SECTION: <i>Oral Solids</i>	PRODUCT PREPARATION	PRODUCT COLLECTION
Date: 11/02/01 Signature: [signature]	Date: 12-04-01 Operator's signature: [signature] Verifier's signature: [signature] Chief's signature: [signature]	Date: 12-6-01 Signature: [signature]
QUALITY ASSURANCE APPROVAL OR DELEGATION FOR THE COLLECTION OF THE PRODUCT STILL UNDER ANALYSIS		
Date: 12/4/2001 Signature: [signature]		

MTH014 4

*\*Delivered all of it to stock*  
[signature] 12 April 2001



LOT 183G02

**PHARMACEUTICAL DEVELOPMENT  
EXCIPIENTS REQUEST**

PRODUCT: <i>SU 10398</i>				
LOT/PREPARATION: <i>183G02</i>				
PHARMACEUTICAL FORM: <i>CPS</i>			DOSAGE: <i>50 mg (AS FREE BASE)</i>	
SCOPE OF THE PREPARATION: <i>CLINICAL SUPPLY</i>				
EXCIPIENT NAME	CODE	LOT	QUANTITY (in grams)	UNDER ANALYSES
<i>F3 T/C SWEDISH ORANGE GEL CAPSULES</i>	<i>1491</i>	<i>AE283</i>	<i>CIRCA 65000 CPS</i>	<i>*</i>
			<i>= p<sub>0</sub>2.3185</i>	
		<i>[initials] 03/04/01</i>		
NOTES: <i>MAKE READY BEFORE 06 APRIL 2001 [initials]</i>				
REQUESTING SECTION: <i>ORAL SOLIDS</i>		PRODUCT PREPARATION		PRODUCT COLLECTION
Date: <i>03/04/01</i>		Date: <i>04-04-01</i>		Date: <i>12-4-01</i>
Signature: <i>[signature]</i>		Operator's signature: <i>[signature]</i>		Signature: <i>[signature]</i>
		Verifier's signature: <i>[signature]</i>		
		Chief's signature: <i>[signature]</i>		

MTF017\_5

WAREHOUSING Section

**PRODUCT**

**F3 T/C SWEDISH ORANGE GEL CPS**

**LOT**

**AE283**

Scale	1
Gross	3.447 kg
Tare	0.262 kg
<b>NET</b>	<b>3.185 kg</b>
Date	04.04.01
Time	11.21.15
Operator	<i>[signature] CANESSA/ANTONINI</i>

**PHARMACEUTICAL DEVELOPMENT  
PACKAGING MATERIALS REQUEST**

PACKAGING MATERIALS					
MATERIAL	CODE	LOT	QUANTITY REQUESTED	QUANTITY SENT	UNDER ANALYSES
KRAFT BARRELS	771350000	VARIOUS	Nº2	2	
PE BAGS FOR BARRELS	735573000	AA39N054	Nº10	10	
PE BAGS 350x580mm	735190000	AA38L198	Nº10	10	
PE BAGS 280x380mm	735170000	AA39D091	Nº15	15	

PRODUCT TO BE PACKAGED: <u>SU10398</u>	PHARMACEUTICAL FORM: <u>CPS</u>
Dosage: <u>50mg (AS FREE BASE)</u>	Lot: <u>183G02.</u>

REQUESTING SECTION <u>ORAL SOLIDS</u>	PRODUCT PREPARATION	PRODUCT COLLECTION
Date: <u>03/04/01</u> Signature: <u>[signature]</u>	Date: <u>05-04-01</u> Operator's signature: <u>[signature]</u> Verifier's signature: Chief's signature: <u>[signature]</u>	Date: <u>11-4-01</u> Signature: <u>[signature]</u>

MTH014\_4

NOTE MAKE READY BEFORE 06 APRIL 2001 [initials]

## DETERGENTS/DISINFECTANTS REQUEST FOR THE CLEANING OF THE PLANTS

PRODUCT: SU10398				
LOT/PREPARATION: I83G02				
DETERGENT/DISINFECTANT NAME	CODE	LOT	QUANTITY	UNDER ANALYSES
PYRONEG	I453	AE259	5 Ky	*
NOTES: _____ _____				
REQUESTING SECTION: <u>ORAL SOLIDS</u>	PRODUCT PREPARATION		PRODUCT COLLECTION	
Date: 03/04/01 Signature: [signature] _____	Date: 04-04-01 Operator's signature: [signature] Verifier's signature: [signature] Chief's signature: [signature]		Date: 22-4-01 Signature: [signature]	

MTF017 5

[signature]

## WAREHOUSING Section

PRODUCT

PYRONEG

LOT

AE259

Scale

---

1

Gross

5.263 kg

Tare

0.263 kg

NET

5.000 k

Date

04 04 01

Time

12.35.41

## Opera

12.55.41  
CANESS

LOT 183G02  
ATTACHMENT 8 [initials] 27 MAY 2001



Pharmacia & Upjohn

PILOT PLANT FORMULATION DEVELOPMENT

FINISHED PRODUCT DELIVERY FORM

DATE: 10 / MAY / 01

PRODUCT: SU10398 PREPARATION DATE: 04 / 01 APPROVED  
LOT: 183G02 UNDER ANALYSES X  
DOSAGE: 50 mg FORMULA NO.: 83HCO1

RAW MATERIAL LOT: (A) 5975-MTM-0002-M2 (97.93%)

QUANTITY 41875 + COUNTER SAMPLE 100 TOTAL 41975

ADMINISTRATION: oral ☒ injectable ☐ topical ☐ drops ☐

PHARMACEUTICAL FORM

LYOPHILE ampoule ☐ vial ☐  
SOLUTION/SUSPENSION bottle ☐ vial ☐ ampoule ☐ small flask ☐ bag ☐  
OINTMENT tube ☐ jar ☐  
gel ☐ cream ☐ paste ☐ salve ☐  
TABLET simple ☐ film-coated ☐ sugar-coated ☐  
gastrointestinal ☐ soluble/effervescent ☐

dimensions/form: \_\_\_\_\_

average weight: \_\_\_\_\_

packaging: \_\_\_\_\_

CAPSULE hard gelatin ☒ soft gelatin ☐ 136.1 [initials] 10/5/01

format: 3 average weight: 136.26 µg

color: T/C OPAQUE SWEDISH ORANGE

printing: \_\_\_\_\_

packaging: 1 BARREL

POWDER/GRANULATE oral ☐ injectable ☐ inhalational ☐

packaging: \_\_\_\_\_

STORAGE room temperature ☒ +1°C ☐ -20°C ☐ -80°C ☐

other conditions: \_\_\_\_\_

POSSIBLE NOTES: \_\_\_\_\_

PERSON IN CHARGE: [signature]



**Pharmacia & Upjohn**  
Pharmaceutical Development

**Analyses Report**

Laboratory  
Information  
Management  
System

Submission Id: 115751  
Sample ID (Request No.): 20011710  
User sample ID (Lot): J63G02  
Sample Type: SU010598 (PNU290940A0) CAPSULE 50  
MG

Stability Study: 183G02  
Pack Type: BULK  
Storage Condition: TIME ZERO  
Time Label: TIME ZERO  
Preparation purpose: STABILITY

Condition:  
**APPROVED**  
Status:  
**COMPLETE**

Samples Notes:

Method	Controls	Results	M.U.	Spec?	Lower lim	Upper lim	Book-Data Analyses	Analyst	Verifier
083GCC01#0 1	APPEARANCE	COMPLIANT	bt				75/35 - JUN 2001	PASTO	MARIANI
083GIZ01#01	HPLC IDENTIFICATION	POSITIVE	bt				75/35 - JUN 2001	PASTO	MARIANI
083GT01#01 R	(HPLC TITER) SAMPLE#1	98.3	%		95.000	105.000	75/35 - JUN 2001	PASTO	MARIANI
	(HPLC TITER) SAMPLE#2	97.7	%		95.000	105.000	75/35 - JUN 2001	PASTO	MARIANI
	(HPLC TITER) SAMPLE#3	98.1	%		95.000	105.000	75/35 - JUN 2001	PASTO	MARIANI
	AVERAGE TITER (HPLC)	98.03	%		95.000	105.000	75/35 - JUN 2001	PASTO	MARIANI
083GSC01#00	TOTAL CORRELATED SUBSTANCES	1.3	%				75/35 - JUN 2001	PASTO	MARIANI
	IMPURITIES NOT LESS THAN 0.05%	SEE TEXT	%				75/35 - JUN 2001	PASTO	MARIANI
083GKF01#00	(HUMIDITY) SAMPLE #1	1.48	%				75/35 - JUN 2001	PASTO	MARIANI
	(HUMIDITY) SAMPLE #2	1.63	%				75/35 - JUN 2001	PASTO	MARIANI
	(HUMIDITY) SAMPLE #3	1.67	%				75/35 - JUN 2001	PASTO	MARIANI
	AVERAGE HUMIDITY	1.59	%				75/35 - JUN 2001	PASTO	MARIANI



**Pharmacia & Upjohn**  
Pharmaceutical Development

**Analyses Report**

Laboratory  
Information  
Management  
System

Submission ID: 113751  
Sample ID (Request No.): 20011710  
User sample ID (lot): 183G02  
Sample Type: SU010395 (PNU290840AD) CAPSULE 50  
MG

Stability Study: 183G02  
Pack Type: BULK  
Storage Condition: TIME ZERO  
Time Label: TIME ZERO  
Preparation purpose: STABILITY

Condition:  
**APPROVED**  
Status:  
**COMPLETE**

Samples Notes:

Method	Controls	Results	M.U.	Spec?	Lower lim	Upper lim	Book-Data Analyses	Analyst	Verifier
083GUP01#01	SAMPLE#1	97.4	%		85.000	- 115.000	75/35 - JUN 2001	PASTO	MARIANI
	SAMPLE#2	101.2	%		85.000	- 115.000	75/35 - JUN 2001	PASTO	MARIANI
	SAMPLE#3	101.3	%		85.000	- 115.000	75/35 - JUN 2001	PASTO	MARIANI
	SAMPLE#4	100.3	%		85.000	- 115.000	75/35 - JUN 2001	PASTO	MARIANI
	SAMPLE#5	102.2	%		85.000	- 115.000	75/35 - JUN 2001	PASTO	MARIANI
	SAMPLE#6	98.2	%		85.000	- 115.000	75/35 - JUN 2001	PASTO	MARIANI
	SAMPLE#7	96.4	%		85.000	- 115.000	75/35 - JUN 2001	PASTO	MARIANI
	SAMPLE#8	97.2	%		85.000	- 115.000	75/35 - JUN 2001	PASTO	MARIANI
	SAMPLE#9	95.5	%		85.000	- 115.000	75/35 - JUN 2001	PASTO	MARIANI
	SAMPLE#10	105.6	%		85.000	- 115.000	75/35 - JUN 2001	PASTO	MARIANI
	SAMPLE#11	97.3	%		85.000	- 115.000	75/35 - JUN 2001	PASTO	MARIANI
	SAMPLE#12	97.3	%		85.000	- 115.000	75/35 - JUN 2001	PASTO	MARIANI
	SAMPLE#13	98.2	%		85.000	- 115.000	75/35 - JUN 2001	PASTO	MARIANI
	SAMPLE#14	103.7	%		85.000	- 115.000	75/35 - JUN 2001	PASTO	MARIANI
	SAMPLE#15	99.2	%		85.000	- 115.000	75/35 - JUN 2001	PASTO	MARIANI
	SAMPLE#16	95.7	%		85.000	- 115.000	75/35 - JUN 2001	PASTO	MARIANI
	SAMPLE#17	95.8	%		85.000	- 115.000	75/35 - JUN 2001	PASTO	MARIANI
	SAMPLE#18	99.8	%		85.000	- 115.000	75/35 - JUN 2001	PASTO	MARIANI
	SAMPLE#19	96.3	%		85.000	- 115.000	75/35 - JUN 2001	PASTO	MARIANI



**Pharmacia & Upjohn**  
Pharmaceutical Development

**Analyses Report**

Laboratory  
Information  
Management  
System

Submission Id: 113751  
Sample ID (Request No.): 20011710  
User sample ID (lot): **103G02**

Sample Type: SU010398 (PNU290940AD) CAPSULE 50  
MG

Stability Study: 183G02  
Pack Type: BULK  
Storage Condition: TIME ZERO  
Time Label: TIME ZERO  
Preparation purpose: STABILITY

Condition:  
**APPROVED**  
Status:  
**COMPLETE**

Samples Notes:

Method	Controls	Results	M.U.	Spec?	Lower lim	Upper lim	Book-Data Analyses	Analyst	Verifier
083GUP01#01	SAMPLE#20	97.5	%		85.000	- 115.000	75/35 - JUN 2001	PASTO	MARIANI
	AVERAGE WEIGHT	98.61	%		85.000	- 116.000	75/35 - JUN 2001	PASTO	MARIANI
	STANDARD DEVIATION	2.922	std				75/35 - JUN 2001	PASTO	MARIANI
	C.V.	2.98	%				75/35 - JUN 2001	PASTO	MARIANI
	MIN. DISSOLUTION VALUE	95.5	%		85.000	- 115.000	75/35 - JUN 2001	PASTO	MARIANI
	MAX. DISSOLUTION VALUE	105.6	%		85.000	- 115.000	75/35 - JUN 2001	PASTO	MARIANI
083GDS01#01_12									
083GMC01#01	TOTAL AEROBIC BACTERIA	IN LIMITS	bd				BIOLAB MAY 2001	PASTO	MARIANI
	ESCHERICHIA COLI	ABSENT	bd				BIOLAB MAY 2001	PASTO	MARIANI
	SALMONELLA	ABSENT	bd				BIOLAB MAY 2001	PASTO	MARIANI
	FUNGI	IN LIMITS	bd				BIOLAB MAY 2001	PASTO	MARIANI
	STAPHYLOCOCCUS AUREUS	ABSENT	bd				BIOLAB MAY 2001	PASTO	MARIANI
	ENTEROBACTERIA	ABSENT	bd				BIOLAB MAY 2001	PASTO	MARIANI



**Pharmacia & Upjohn**  
Pharmaceutical Development

**Analyses Report**

Laboratory  
Information  
Management  
System

Submission Id: 113751  
Sample ID (Request No.): 20011710  
User sample ID (lot): J83G02  
Sample Type: SU010328 (PNU290940AD) CAPSULE 50  
MG

Stability Study: I63G02  
Pack Type: BULK  
Storage Condition: TIME ZERO  
Time Label: TIME ZERO  
Preparation purpose: STABILITY

Condition :  
**APPROVED**  
Status :  
**COMPLETE**

Samples Notes:

Method	Controls	Results	M.U.	Spec?	Lower lim	Upper lim	Book-Data Analyses	Analyst	Verifier
083GDS01#01 _12	SAMPLE#1	97.5	%		80.000		75/35 - JUN 2001	PASTO	MARIANI
	SAMPLE#2	100.5	%		80.000		75/35 - JUN 2001	PASTO	MARIANI
	SAMPLE#3	99.0	%		80.000		75/35 - JUN 2001	PASTO	MARIANI
	SAMPLE#4	95.5	%		80.000		75/35 - JUN 2001	PASTO	MARIANI
	SAMPLE#5	99.9	%		80.000		75/35 - JUN 2001	PASTO	MARIANI
	SAMPLE#6	97.2	%		80.000		75/35 - JUN 2001	PASTO	MARIANI
	SAMPLE#7	98.7	%		80.000		75/35 - JUN 2001	PASTO	MARIANI
	SAMPLE#8	99.3	%		80.000		75/35 - JUN 2001	PASTO	MARIANI
	SAMPLE#9	97.8	%		80.000		75/35 - JUN 2001	PASTO	MARIANI
	SAMPLE#10	98.7	%		80.000		75/35 - JUN 2001	PASTO	MARIANI
	SAMPLE#11	99.0	%		80.000		75/35 - JUN 2001	PASTO	MARIANI
	SAMPLE#12	99.0	%		80.000		75/35 - JUN 2001	PASTO	MARIANI
	DISSOLUTION	98.51	%		80.000		75/35 - JUN 2001	PASTO	MARIANI
	STANDARD DEVIATION	1.335	std-1				75/35 - JUN 2001	PASTO	MARIANI
	C.V.	1.36	%				75/35 - JUN 2001	PASTO	MARIANI
	MIN. DISSOLUTION VALUE	95.5	%		80.000		75/35 - JUN 2001	PASTO	MARIANI
	MAX. DISSOLUTION VALUE	100.5	%		80.000		75/35 - JUN 2001	PASTO	MARIANI





**Pharmacia & Upjohn**  
Pharmaceutical Development

**Analyses Report**

Laboratory  
Information  
Management  
System

Submission Id: 113751  
Sample ID (Request No.): 20011710  
User sample ID (Lot): **183G02**  
Sample Type: SUD10396 (PNU290940AD) CAPSULE 50  
MG

Stability Study: 183G02  
Pack Type: BULK  
Storage Conditions: TIME ZERO  
Time Label: TIME ZERO  
Preparation purpose: STABILITY

Condition:  
**APPROVED**  
Status:  
**COMPLETE**

**Samples Notes:**

Observations on the Method 063GCO1#00, component: APPEARANCE  
Opaque brick red colored shaped hard gelatin capsule containing yellow orange colored powder

Observations on the Method 063GDS01#01\_12, component: MAXIMUM DISSOLUTION VALUE  
no comment

Observations on the Method 063GDS01#01\_12, component: DISSOLUTION  
Average dissolution at 15 minutes: 97.7% (C.V. 2.6%)  
Average dissolution at 30 minutes: 96.5% (C.V. 1.4%)  
Average dissolution at 45 minutes: 99.5% (C.V. 1.7%)

Observations on the Method 063GMC01#01, component: TOTAL AEROBIC BACTERIA  
Total bacteria count < 10 Ufc/g (Report No. 01.1231)

Observations on the Method 063GMC01#01, component: ENTEROBACTERIA  
Enterobacteria < 10 Ufc/g (Report No. 01.1231)

Observations on the Method 063GMC01#01, component: ESCHERICHIA COLI  
Test done on 10g (Report No. 01.1231)

Observations on the Method 063GMC01#01, component: FUNGI  
Fungi (yeasts and molds) < 10 Ufc/g (Report No. 01.1231)

Observations on the Method 063GMC01#01, component: SALMONELLA  
Test done on 10g (Report No. 01.1231)

Observations on the Method 063GMC01#01, component: STAPHYLOCOCCUS AUREUS  
Report No. 01.1231

Observations on the Method 063GSC01#00, component: IMPURITIES NOT LESS  
THAN 0.05%



**Pharmacia & Upjohn**  
Pharmaceutical Development

**Analyses Report**

Laboratory  
Information  
Management  
System

Submission Id: 113751  
Sample ID (Request No.): 20011710  
User sample ID (lot): **I83G02**  
Sample Type: SU010398 (PNU290940AD) CAPSULE 50  
MG

Stability Study: I83G02  
Pack Type: BULK  
Storage Condition: TIME ZERO  
Time Label: TIME ZERO  
Preparation purpose: STABILITY

Condition :  
**APPROVED**  
Status :  
**COMPLETE**

**Samples Notes:**

Observations on the Method 063GSC01#00, component: IMPURITIES NOT LESS  
THAN 0.05%

RRT 0.690 (Impurities A): 0.16%  
RRT 0.913 (Impurities B): 0.09%  
RRT 0.926 (Impurities C): 0.15%  
RRT 1.069 (Impurities D): 0.28%  
RRT 1.089 (Impurities E): 0.11%  
RRT 1.142 (Impurities F): 0.07%  
RRT 1.261 (Impurities G): 0.13%  
RRT 1.272 (Impurities H): 0.15%  
RRT 1.646 (Impurities I): 0.16%

Observations on the Method 063GUP01#00, component: AVERAGE WEIGHT  
The average weight of 98.51% corresponds to 87.65 mg/capsule

Approved by: MARIANI  
Approval date: 13-Jun-2001

<b>PRODUCT/MATERIAL:</b> <b>SU010398 50mg CAPSULE</b>		<b>NUMBER(S) OF THE LOT(S)/PACKAGING RUN:</b>  <b>I83K01 and I 83G02</b>	
<b>DATE (AND TIME WHERE POSSIBLE) OF THE DEVIATION:</b> <b>15 June 2001</b>			
<b>DESCRIPTION OF PROPOSED DEVIATION:</b>			
1 - Absence of the following instructions on the Processing Sheet regarding what is described in the product's Process Flow Chart: start time and end time registration of lot I82K01, Granulation operations (start time Num. 3) and Num. 13 and on the I82G02 sheet, operations Num. 2 and Num. 12. 2 - Encapsulation operations during the period from the 12th to the 24th of April			
		<b>SIGNATURE/DATE:</b> [signature] June 15th, 2001	
<b>REASON FOR THE DEVIATION AND INVESTIGATION:</b>			
1 - lack of registration indications in the filling in the Processing Sheet by Oral Solids 2 - The extended period is due to the Easter Holidays and unexpected technical difficulties, namely the adhesion of the mix to the walls of the hopper. Stability data for the month regarding a technical capsule lot (1985-013) support the granulate's stability for a period of 12 days, which is the duration of the capsule preparation process.			
		<b>SIGNATURE/DATE:</b> [signature] June 15th, 2001	
<b>IMMEDIATE CORRECTIVE ACTIONS:</b> <b>NONE</b>			
<b>SIGNATURE/DATE:</b>			
<b>SECTION CHIEF:</b>			
<b>QA/QUALITY SYSTEMS:</b>			
<b>1 - Oral Solids will perform the filling out of the Process Sheet in accordance with the Process Flow Chart prior to the next productions</b> <b>2 - If in future production the time is prolonged longer than expected, appropriate analytic controls will be performed on the granulate to verify the product's stability. For this purpose a portion of granulate must be collected at the beginning of encapsulation, which is to be kept at +2° - +8° C, in case an initial reference sample is needed.</b>			
		<b>SIGNATURE/DATE:</b> [signature] June 15th, 2001	
<b>SECTION CHIEF:</b>			
<b>QA/QUALITY SYSTEMS:</b>			
[signature]		[signature] 18/06/01	

COMPLAINT N° 24 / 2002

Product: SU011248 L-malate 50 mg capsules - Lot 183C02 (manufacturing date April 2001)

REASON FOR THE COMPLAINT:

EVIDENCED IN A CLINICAL CENTER A CAPSULE WITH CRUSHED END CONTAINED IN ONE BOTTLE OF THE DRUG  
(COMMUNICATION FROM SUGEN - see enclosure 1)

ENCHARGED PERSON OF THE INVOLVED SECTION

ALESSANDRA CAVALLO

DATE OF RECEIPT OF THE COMPLAINT

- ORAL SOLIDS R&D: 21.11.2002
- QA: 21.11.2002

DETAILED DESCRIPTION OF THE DEFECT:

PRIMARY DEFECT: CAPSULE WITH CRUSHED END

INVESTIGATION RESULTS:

THE 3% OF THE LOT WAS INITIALLY INSPECTED (1260 CPS ON 41985), FINDING OUT 7 UNITS WITH BROKEN END. AS THE ACCEPTANCE LIMIT IS OF THREE UNITS, THE LOT WAS 100% VISUALLY INSPECTED (MANUAL OPERATION) ACCORDING TO SOP SF.TF 025, DISCHARGING ALL THE DEFECTIVE CAPSULES (6 UNITS FOR BROKEN END, 3 UNITS FOR BROKEN BODY, 3 UNITS FOR VISUALIZATION OF THE BODY ON THE TOP AND 378 UNITS FOR CRUSHED ENDS)  
THE ROOT CAUSE OF THE EVIDENCED DEFECTS IS ATTRIBUTABLE TO THE NOT CORRECT REGULATION OF THE PUNCHES CLOSING THE CAPSULES

IMMEDIATE CORRECTIVE ACTIONS (indicate who is responsible for / completion date):

NO IMMEDIATE ACTION - THE EVIDENCED DEFECT DO NOT HAVE ANY IMPACT ON THE PRODUCT EFFICACY OR SAFETY FOR THE PATIENT

SIGNATURE/DATE: *Luigi Facchetti* 18/12/2002

MEAN/LONG TERM CORRECTIVE ACTIONS (indicate who is responsible for / completion date):

- THE PUNCHES CLOSING THE CAPSULES ARE SET UP BEFORE EACH PRODUCTION, TRYING TO OPTIMIZE THE REGULATION. PARTICULAR ATTENTION WILL BE DEVOTED TO THIS ASPECT AS FOR FURTHER ENCAPSULATION PROCESSES.  
(Alessandra Cavallo Responsible for).
- A TRAINING SESSION FOR THE OPERATORS PERFORMING VISUAL INSPECTIONS WILL BE ORGANIZED  
(Alessandra Cavallo Responsible for; foreseen completion date: end February, 2003).

CONCLUSIONS:

THE EVIDENCED DEFECT IS ATTRIBUTABLE TO OPERATOR'S OVERLOOKING. NO COMPLAINT AS FOR THE SAME DEFECT HAS UNTIL NOW BEEN RECEIVED.

SIGNATURE/DATE: *Luigi Facchetti*

18/12/2002

**REASON FOR THE COMPLAINT:**

**A CLINIC REPORTED THE PRESENCE OF A CAPSULE WITH A CRUSHED BASE IN A PACKAGE OF THE DRUG (SUGEN COMMUNICATION- ATTACHMENT 1)**

**ENCHARGED PERSON OF THE INVOLVED SECTION**

**ALESSANDRA CAVALO**

**CLAIM RECEIVED ON AND BY:**

- 21 November 2002 by ORAL SOLIDS R&D
- 21 November 2002 by QA

**DETAILED DESCRIPTION OF THE DEFECT:**

**PRIMARY DEFECT – CAPSULE BASE CRUSHED –**

**INVESTIGATION RESULTS:**

**3% OF THE LOT WAS INITIALLY INSPECTED (1260 CAPSULES OUT OF 41985) WITH THE DETECTION OF 7 DEFECTIVE UNITS DUE TO RUPTURE AT THE TIP. SINCE THE LIMIT OF THREE DEFECTIVE UNITS WAS EXCEEDED, THE LOT WAS SORTED BY UNIT (MANUAL OPERATION), ACCORDING TO THE PROCEDURE SF.TF 025, IN ORDER TO ELIMINATE ALL DEFECTIVE CAPSULES (6 UNITS FOR TIP RUPTURE AT THE TIPS, 3 UNITS FOR FRACTURE OF THE CAPSULE BODY, 3 UNITS FOR THE VISUALIZATION OF THE BODY ON THE TOP, AND 378 UNITS FOR CRUSHED TIPS). THE CAUSE OF THE DEFECTS FOUND CAN BE ATTRIBUTED TO THE INCORRECT REGULATION OF THE PUSHING ELEMENTS FOR CLOSING THE CAPSULES.**

**IMMEDIATE CORRECTIVE ACTIONS (indicate who is responsible for / completion date):**

**NONE- The nature of the discovered defects did not have any impact on patient safety and/or the drugs efficacy.**

**SIGNATURE/DATE:** [signature] 18/12/2002

**MEAN/LONG TERM CORRECTIVE ACTIONS (indicate who is responsible for / completion date):**

**- THE PUSHING ELEMENTS TO CLOSE THE CAPSULES ARE MOUNTED PRIOR TO EACH ENCAPSULATION, IN AN ATTEMPT TO ACHIEVE THE OPTIMAL REGULATION. PARTICULAR ATTENTION WILL BE GIVEN AT THE START OF FUTURE PROCESSING.  
(Person in Charge Alessandra Cavallo).**

**- TRAINING SESSIONS WILL BE ARRANGED FOR THE OPERATORS OF THE SORTING OPERATIONS  
(Person in Charge Alessandra Cavallo; date planned for completion: End of February 2003).**

**CONCLUSIONS:**

**THE DETECTED DEFECT IS ATTRIBUTED TO THE OPERATOR'S MISTAKE DURING THE SORTING PHASE. AS YET NO CLAIM OF ANY KIND HAS BEEN RECEIVED.**

**SIGNATURE/DATE:** [signature]  
18/12/2002